

Instructions for Use: Humipak

Brand Name of Product	Humipak
Generic Name of Product	Transport pouch
Product Code Number(s)	HPSS2035-25, HPSS2035-100, HPSS4243-25, HPSS4243-100, HPSS4267-25, HPSS4267-100, HPSS6577-10, HPSS6577-50, HPSS2035, HPSS4243, HPSS4267, HPSS6567, HPSS6577
Purpose of Product	To provide a humid environment delaying the drying of contaminants during transport Prior to any required cleaning, disinfection, reprocessing and/or sterilization.
Range of Applications for Product	Contaminated items requiring further processing.
Key Specifications of Product	<ul style="list-style-type: none"> • Self-seal strip. • A special absorbent layer which when wetted with water creates a humid atmosphere within the pouch. • Two layers of transparent film (containing the special absorbent layers) form the pouch. • The transparent film (in addition to maintaining the humid atmosphere) allows the contents to be clearly observed. <ul style="list-style-type: none"> ○ Four (4) sizes: a) Small, b) Half-Din, c) Full-Din, d) Double-Din • The Humipak is NOT intended for: <ul style="list-style-type: none"> ○ <i>Patient use</i> (e.g., not to be used in the diagnosis of disease or other conditions, or in the cure, mitigation, treatment, or prevention of disease). ○ <i>Providing any element for cleaning, disinfection, or sterilization of reusable medical devices</i>, which require cleaning, disinfection, or sterilization prior to reuse. ○ <i>Sterile barrier packaging</i> for medical devices.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Store in a clean, dry place at room temperature.
Packaging Contents	N/A
Shelf Life	Five (5)-year shelf life from the date of manufacture.

Instructions for Using Product	
Description of Use(s)	Placing items into the Humipak after use will provide a humid environment delaying the drying of contaminants during transport Prior to any required cleaning, disinfection, reprocessing and/or sterilization.
Preparation for Use	Consult the manufacturer's IFU for pre-treatments at the point of use.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Pick the correct size of Humipak pouch for the items that are to be kept moist until reprocessing. 2. Add the specified amount of water per the directions found on the Humipak pouch. 3. Water will wick and be absorbed throughout the Humipak pouch because of the special absorbent layer between the transparent films. 4. Place items in the pouch. 5. Peel off the protective strip from the adhesive band. 6. Fold over the adhesive strip and press down firmly from the center outwards to ensure a watertight seal. 7. Send sealed pak to the appropriate area for reprocessing to begin as soon as possible. 8. Cut open—preferably cutting along the side or end farthest away from the body to avoid direct contact with air expelled from the pouch. 9. Ensure environmental conditions within the decontamination area are controlled to minimize the spread of airborne organisms.

Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	<ul style="list-style-type: none"> • Items placed in the Humipak should be reprocessed as soon as possible. • Single-use only. • Do not reuse or sterilize the Humipak. • Use proper PPE when opening the Humipak. • The Humipak is NOT a sharps container and is NOT puncture proof. • Follow your facility's policy for handling contaminated reusable sharps in compliance with OSHA Guidelines.
Disposal	Treat used Humipak as a biohazard and dispose of in accordance with facility's policy for contaminated items.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	SST
Other Product Support Documents	SST Instrument Retrieval Brochure, SST Instrument Retrieval Price List, Humipak Pocket Brochure
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com