

Brand Name of Product	Humipak
Generic Name of Product	Transport pouch
Product Code Number(s)	HPSS2035-25, HPSS2035-100, HPSS4243-25, HPSS4243-100, HPSS4267-25, HPSS4267-100, HPSS6577-10, HPSS6577-50, HPSSR01-100, HPSSR01-25, HPSSR04-100, HPSSR04-25
Purpose of Product	Provides a humid environment delaying the drying of contaminants during transport prior to required cleaning, disinfection, reprocessing and/or sterilization.
Range of Applications for Product	<p>Holds contaminated items and maintains moisture to stop clinical soils from drying before further processing.</p> <p>The Humipak is NOT intended for:</p> <ul style="list-style-type: none"> ○ <i>Patient use</i> (e.g., not to be used in the diagnosis of disease or other conditions, or in the cure, mitigation, treatment, or prevention of disease). ○ <i>Providing any element for cleaning, disinfection, or sterilization of reusable medical devices</i>, which require cleaning, disinfection, or sterilization prior to reuse. ○ <i>Sterile barrier packaging</i> for medical devices.
Key Specifications of Product	<ul style="list-style-type: none"> • Self-seal strip • A special absorbent layer (when wet with water) creates a humid atmosphere within the pouch. • Two (2) layers of transparent film (containing the special absorbent layers) form the pouch. • The transparent film (in addition to maintaining the humid atmosphere) allows the contents to be clearly observed. • Six (6) sizes: <ol style="list-style-type: none"> a. (Item #HPSS2035) Small 200- x 350 mm (7.87- x 13.78 in) b. (Item #HPSS4243) Half-Din 420- x 430mm (16.54- x 16.93 in) c. (Item # HPSS4267) Full-Din 420- x 670 mm (16.54- x 26.38 in) d. (Item # HPSS6577) Double-Din 650- x 770 mm (25.59- x 30.31 in) e. Robotic instruments One and Four Arm Configurations (Item codes): <ul style="list-style-type: none"> ■ HPSSR01-100 ■ HPSSR01-25 ■ HPSSR04-100 ■ HPSSR04-25

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Store in a clean, dry place at room temperature (i.e., 20- to 25 °C [68- to 77 °F]).
Packaging Contents	N/A
Shelf Life	Five (5)-year shelf life from the date of manufacture.

Instructions for Using Product	
Description of Use(s)	Placing used items into the Humipak provides a humid environment and will delay drying of contaminants and clinical soil during transport prior to any required cleaning, disinfection, reprocessing, and/or sterilization.
Preparation for Use	Consult the manufacturer's IFU for pretreatments at the point of use.
Diagrams (drawings, pictures)	




Figure 1



Figure 2



Figure 2a



Figure 3



Figure 4




Figure 5









Figure 5A

		
	Figure 6	Figure 6a
		
	Figure 6b	
		
	Figure 7	Figure 7b
Steps for Use of Product	<ol style="list-style-type: none"> 1. Pick the <i>correct size</i> of Humipak pouch for items to be kept moist until reprocessing. (Fig. 1). 2. Add the <i>specified amount of water</i> per the directions found on the Humipak pouch. (Fig. 2, 2a). 3. Water will wick and be absorbed throughout the Humipak pouch because of the special absorbent layer between the transparent films. (Fig. 3). 4. Place item(s) in the pouch. (Fig.4). 5. Peel off the protective strip from the adhesive band. (Fig. 5, 5a). 6. Fold over the adhesive strip and press down firmly from the center outwards to ensure a watertight seal. (Fig. 6, 6a, 6b). 7. Send sealed Humipak to the appropriate area for reprocessing to begin as soon as possible. 8. <i>Cut open</i>—preferably cutting along the side or end farthest away from the body to avoid direct contact with air expelled from the pouch. (Fig. 7, 7b). 9. <i>Ensure controlled environmental conditions</i> within the decontamination area to minimize the spread of airborne organisms. 	
Interpretation of Results	N/A	
Contraindications of Test Results	N/A	
Documentation	N/A	
Special Warnings and Cautions	<ul style="list-style-type: none"> • Items placed in the Humipak should be reprocessed as soon as possible. • Single use only. • Humipak is NOT: <ul style="list-style-type: none"> ○ A sharps container. ○ Puncture proof. • ⚠ Do NOT reuse or sterilize the Humipak. • ⚠ Use proper PPE when opening the Humipak. 	

	<ul style="list-style-type: none"> • ☒ Treat used Humipak(s) as biohazardous. • ⚠ Follow your facility's policy for handling contaminated reusable sharps in compliance with OSHA Guidelines.
Disposal	Treat used Humipak as a ☒ biohazard and dispose of in accordance with facility's policy for contaminated items.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	SST
Other Product Support Documents	SST Instrument Retrieval Brochure, Humipak Pocket Brochure
Reference Documents	N/A
Customer Service Contact	Healthmark, A Getinge company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com