



Instructions for Use: Sample Extraction Tools

Brand Name of Product	Sample Extraction Tool
Generic Name of Product	Single use sample extraction brushes
Product Code Number(s)	CKB-200-200, CKB-230-400, CKB-230-600
Intended Use	The sample extraction tool is intended for extracting samples to test for bacterial growth.
Range of Applications for Product	Extracting samples for testing.
Key Specifications of Product	<ul style="list-style-type: none"> • Proper length and diameters for contact with the internal surfaces. <ul style="list-style-type: none"> - CKB 200-200: 200cm Length x 2mm Diameter - CKB 230-400: 230cm Length x 4mm Diameter - CKB 230-600: 230cm Length x 6mm Diameter

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To extract samples for testing bacterial growth.
Preparation for Use	<ol style="list-style-type: none"> 1. Select the diameter and length for the internal surface/channel. 2. Remove the extraction tool from the packaging.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Insert the sample extraction tool into the channel. 2. Continue to feed the extraction tool through the length of the channel until it reaches the distal end. Never force the extraction tool, to avoid possible damage. 3. Once the extraction tool exits the distal end, cut the entire head of the bristle portion of the extraction tool into the sample collection container. 4. Carefully pull the entire sample extraction tool completely out of the distal end. 5. Dispose of the extraction tool after use.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	<ul style="list-style-type: none"> • Select the extraction tool based upon the diameter and length of the channel chosen. • Do not force the extraction tool into a channel to avoid possible damage.
Disposal	Treat used sample extraction tools as a biohazard and dispose of in accordance with facility policy of contaminated items.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A

Additional Information	N/A
Related Healthmark Products	Flexible Endoscope Sampling Kits
Other Product Support Documents	ProFormance™ Brochure, ProFormancew™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com