

## Instructions for Use: Sampling Swab Kit

<b>Brand Name of Product</b>	Sampling Swab Kit												
<b>Generic Name of Product</b>	Sampling Swab Kit												
<b>Product Code Number(s)</b>	SWK-001, SWK-002, SWK-003, SWK-004, SWK-005, SWK-006, SWK-006, SWK-007, SWK-008, SWK-009, SWK-010, SWK-011, SWK-012												
<b>Intended Use</b>	To obtain a sample to be used for residual testing.												
<b>Range of Applications for Product</b>	Obtaining a sample by a swab for residual testing.												
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• 12 swabs</li> <li>• 12 prepackaged vials of water</li> <li>• 1 pair of scissors</li> </ul> <table> <tr> <td>SWK-001: 1.1mm x 305mm</td> <td>SWK-007: 1.1mm x 2000mm</td> </tr> <tr> <td>SWK-002: 2.0mm x 305mm</td> <td>SWK-008: 2.0mm x 2300mm</td> </tr> <tr> <td>SWK-003: 2.7mm x 305mm</td> <td>SWK-009: 2.7mm x 2300mm</td> </tr> <tr> <td>SWK-004: 3.5mm x 305mm</td> <td>SWK-010: 3.5mm x 2300mm</td> </tr> <tr> <td>SWK-005: 4.2mm x 305mm</td> <td>SWK-011: 4.2mm x 2300mm</td> </tr> <tr> <td>SWK-006: 4.7mm x 305mm</td> <td>SWK-012: 4.7mm x 2300mm</td> </tr> </table>	SWK-001: 1.1mm x 305mm	SWK-007: 1.1mm x 2000mm	SWK-002: 2.0mm x 305mm	SWK-008: 2.0mm x 2300mm	SWK-003: 2.7mm x 305mm	SWK-009: 2.7mm x 2300mm	SWK-004: 3.5mm x 305mm	SWK-010: 3.5mm x 2300mm	SWK-005: 4.2mm x 305mm	SWK-011: 4.2mm x 2300mm	SWK-006: 4.7mm x 305mm	SWK-012: 4.7mm x 2300mm
SWK-001: 1.1mm x 305mm	SWK-007: 1.1mm x 2000mm												
SWK-002: 2.0mm x 305mm	SWK-008: 2.0mm x 2300mm												
SWK-003: 2.7mm x 305mm	SWK-009: 2.7mm x 2300mm												
SWK-004: 3.5mm x 305mm	SWK-010: 3.5mm x 2300mm												
SWK-005: 4.2mm x 305mm	SWK-011: 4.2mm x 2300mm												
SWK-006: 4.7mm x 305mm	SWK-012: 4.7mm x 2300mm												

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	<ul style="list-style-type: none"> <li>• Prepackaged vials of water expire 2 years from date of manufacture.</li> <li>• Expiration date is printed on the package.</li> </ul>

Instructions for Using Product	
<b>Description of Use(s)</b>	Obtain sample after cleaning and prior to further processing to test for residuals.
<b>Preparation for Use</b>	<ul style="list-style-type: none"> <li>• Sampling of the item is done after cleaning prior to further processing.</li> <li>• Select the appropriate size sample swab.</li> <li>• To avoid contamination that might lead to false positive test results, do not touch the swab.</li> </ul>
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Wet surfaces are sampled with a dry swab.</li> <li>2. Dry surfaces are sampled by moistening the swab with a drop of prepackaged, chlorine-free water (supplied).</li> <li>3. Swab the target area of the item.</li> <li>4. For lumen surface sampling, if the design of the item allows it (e.g., open ended), insert the swab end into the item. Push it all the way through one (1) time.</li> <li>5. Cut the swab end of the shaft with scissors (supplied) directly into the reagent vial.</li> <li>6. In the case of open-ended lumen surface testing, remove the remaining segment by pulling the shaft out the distal end.</li> </ol>
<b>Interpretation of Test Results</b>	Consult the IFU of the reagent test for interpretation of results.
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	Dispose of in a biohazard container, according to facility policy.

Reprocessing Instructions	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	<ul style="list-style-type: none"> <li>• The sampling swabs are single-use.</li> </ul>

	<ul style="list-style-type: none"> <li>The scissors may be reused. Clean the scissors in between tests with an IPA wipe.</li> </ul>
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	
<b>Related Healthmark Products</b>	PT-202, HC-101
<b>Other Product Support Documents</b>	ProFormance™ Brochure, ProFormance™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com