

Instructions for Use: Clamps & Spreaders

Brand Name of Product	Clamps & Spreaders
Generic Name of Product	Clamps & Spreaders
Product Code Number(s)	BC-01, BC-02, BC-03, BC-04, SSIS-001, T0200
Purpose of Product	To keep an item in the open position while being cleaned and sterilized.
Range of Applications for Product	<ul style="list-style-type: none"> Clamps keep joints and jaws of the bone forceps open. Spreaders are for single and double action jointed items unable to be dismantled with or without ring handles. The jaws remain open during cleaning and sterilization processes.
Key Specifications of Product	<ul style="list-style-type: none"> Stainless Steel Silicone Covering

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Keeps an item in the open position while being cleaned and sterilized.
Preparation for Use	Select the appropriate size clamp before using.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> Open item to the appropriate width. Insert ends of the clamps or spreaders through the handles' loops and rings or squeeze handles to expose joint and place clamp around both handles to hold joint open. Continue normal procedures according to the manufacturer's guidelines.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	AAMI, AORN, device manufacturers, and other industry experts strongly recommend items should be kept in an open position for maximum exposure.
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	Can be cleaned in an automatic washer with other items.
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> Reprocessing steam sterilization: <ul style="list-style-type: none"> 132 °C (270 °F) for four (4) minutes 135 °C (275 °F) for three (3) minutes. Steam sterilization with instruments: <ul style="list-style-type: none"> 132 °C (270 °F) for four (4) minutes 135 °C (275 °F) for three (3) minutes. Compatible with EO sterilization.
Storage	N/A

Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com