



Brand Name of Product	Container Loks with Indicator Dot
Generic Name of Product	Rigid sterilization container tamper-evident seal with process indicator
Product Code Number(s)	3030 BLWDTS, 3030 YLWDTS, 3030 ORWDTS, 4040 GNWDTS, 5030 RDWDTS, 6020 WTWDTS, 605LNG GRWDTS, 605M GRWDTS, 9393 RDWDTS, 9393 WTWDTS
Purpose of Product	Seal closed sterilization container(s) with a tamper-evident lock to indicate the container(s) have not been opened since sterilization.
Range of Applications for Product	N/A
Key Specifications of Product	Blue, Yellow, Orange, White, Red, Green (4040 GNWDTS only available in green)

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	<ul style="list-style-type: none"> ● Store in a dry condition at room temperature <ul style="list-style-type: none"> ○ (< 60% rh) ○ 15- to 30 °C (59- to 86 °F). ● Do not store near strong alkaline or acidic products (e.g., cleaning/disinfecting agents). ● Store away from chemicals.
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Seal rigid sterilization container with a tamper-evident lock to indicate the container has not been opened since sterilization.
Preparation for Use	N/A
Diagrams (drawings, pictures)	 
Steps for Use of Product	<ol style="list-style-type: none"> 1. Pass arrow shackle through opening of item to be sealed. 2. Firmly push arrow shackle into locking device until it engages. 3. Pull on seal to ensure it has been properly locked.
Interpretation of Results	<i>Steam sterilization:</i> Reactive ink will turn gray/brown/black when exposed to the time/temperature of common steam sterilization cycles (see “Sterilization” section p. 2 below).
Contraindications of Test Results	N/A

Documentation	N/A
Special Warnings and Cautions	If a seal is broken, assume sterility has been compromised and the tray must be reprocessed.
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> • Steam indicator dot (Type 1 chemical indicator), it is suitable for typical hospital steam sterilization cycles, including gravity: <ul style="list-style-type: none"> ○ 121 °C (250 °F) for 30 minutes ○ 132 °C (270 °F) for 4 minutes ○ 135 °C (275 °F) for 3 minutes Note: The indicators will not change in dry heat. • Maximum time/temperature tested: 135 °C for 30 minutes.
Storage	<ul style="list-style-type: none"> • Store in a dry condition at room temperature <ul style="list-style-type: none"> ○ (< 60% rh) ○ 15- to 30 °C (59- to 86 °F). • Do not store near strong alkaline or acidic products (e.g., cleaning/disinfecting agents). • Store away from chemicals.
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Sterilization Packaging Products Brochure, Sterilization Packaging Products Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com