






## Instructions for Use: Crescent Style Stringer

|  |  |
|--|--|
| <b>Brand Name of Product</b>             | Crescent Style Stringers   |
| <b>Generic Name of Product</b>           | Crescent Style Instrument Stringers  |
| <b>Product Code Number(s)</b>            | ST2.64C, ST2.84C, ST2.124C, ST2.144C, ST2.164C, ST2.184C, ST3-64C, ST3-84C, ST3-104C, ST3-124C, ST3-144C, ST3-164C, ST3-184C |
| <b>Purpose of Product</b>                | This product is used to expand instruments to better expose hard-to-reach areas during cleaning and sterilization.           |
| <b>Range of Applications for Product</b> | N/A  |
| <b>Key Specifications of Product</b>     | Manufactured of stainless steel  |

| Shipping & Storage                            |     |
|---|-----|
| <b>Shipping Conditions &amp; Requirements</b> | N/A |
| <b>Storage Conditions</b>                     | N/A |
| <b>Packaging Contents</b>                     | N/A |
| <b>Shelf Life</b>                             | N/A |

| Instructions for Using Product       |   |
|--------------------------------------|---|
| <b>Description of Use(s)</b>         | Keeps instruments maximally exposed for hard-to-reach areas during cleaning and sterilization.  |
| <b>Preparation for Use</b>           | N/A   |
| <b>Diagrams (drawings, pictures)</b> |   <br>    |
| <b>Steps for Use of Product</b>      | <ol style="list-style-type: none"> <li>1. Place one hand at each end of the stringer.</li> <li>2. Using the hand over the crescent-end of the stringer, pull to separate the two pieces.</li> <li>3. String handled instruments (in set assembly order) onto the stringer by placing one ring on each of the stringer arms.             <ol style="list-style-type: none"> <li>a. Continue until all instruments are loaded onto the stringer.</li> <li>b. Lock the stringer by reattaching the crescent-end. (NOTE: You may need to gently pull one of the arms outward to thread into the crescent end.)</li> </ol> </li> </ol> |

|  |     |
|--|-----|
| <b>Interpretation of Results</b>         | N/A |
| <b>Contraindications of Test Results</b> | N/A |
| <b>Documentation</b>                     | N/A |
| <b>Special Warnings and Cautions</b>     | N/A |
| <b>Disposal</b>                          | N/A |

| <b>Reprocessing Instructions</b>            |   |
|---|---|
| <b>Point of Use</b>                         | N/A   |
| <b>Preparation for Decontamination</b>      | N/A   |
| <b>Disassembly Instructions</b>             | N/A   |
| <b>Cleaning – Manual</b>                    | <ul style="list-style-type: none"> <li>Before initial use of the stringer, cleaning is required before initial use of the stringer. <ul style="list-style-type: none"> <li>Manually wash using neutral pH detergent and warm water. (NOTE: Expand and collapse the adjustable assembly area and actuate locking mechanisms during the washing process.)</li> <li>Place the stringer(s) in ultrasonic cleaner for effective cleaning</li> <li>Next place stringer(s) in automated washer system for effective cleaning</li> </ul> </li> <li>After initial use, the stringer is adequately cleaned during the cleaning of the instruments.</li> </ul> |
| <b>Cleaning – Automated</b>                 | <ul style="list-style-type: none"> <li>See instructions under “Steps for Use of Product” for how to string instruments.</li> <li>Consult the instrument manufacturer’s (Mfr.’s) IFU for specific instructions.</li> </ul>   |
| <b>Disinfection</b>                         | N/A   |
| <b>Drying</b>                               | N/A   |
| <b>Maintenance, Inspection, and Testing</b> | Apply instrument milk to sliding crossbar mechanism as needed to ensure smooth action.  |
| <b>Reassembly Instructions</b>              | N/A   |
| <b>Packaging</b>                            | N/A   |
| <b>Sterilization</b>                        | <ol style="list-style-type: none"> <li>See instructions under “Steps for Use of Product” for directions on how to string instruments.</li> <li>When using for sterilization of instruments: <ol style="list-style-type: none"> <li>Open to the optimal width for use in the instrument set</li> <li>Consult the instrument Mfr.’s IFU for specific instructions.</li> </ol> </li> <li>Stringers are materially compatible with steam and EO sterilization.</li> </ol>   |
| <b>Storage</b>                              | N/A   |
| <b>Additional Information</b>               | N/A   |

|  |   |
|--|---|
| <b>Related Healthmark Products</b>     | Stringers   |
| <b>Other Product Support Documents</b> | ProSys™ Brochure, ProSys™ Price List  |
| <b>Reference Documents</b>             | N/A   |
| <b>Customer Service Contact</b>        | Healthmark Industries Company, Inc.<br>18600 Malyn Blvd.<br>Fraser, MI 48026<br>1-586-774-7600<br><a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a><br>hmark.com |