








**Instructions for Use: Magic Touch Magnifier**

<b>Brand Name of Product</b>	Magic Touch Magnifier
<b>Generic Name of Product</b>	Magnifier
<b>Product Code Number(s)</b>	RNDMAG-1, RNDMAGSTND-1
<b>Purpose of Product</b>	A must-have tool with a small size that makes it convenient to carry at all times. Designed to allow for easy and thorough inspection of items.
<b>Range of Applications for Product</b>	Items that need up close visual inspection for possible damage or residual soiling after reprocessing.
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• 3x magnification</li> <li>• Lightweight .45 lbs and durable design</li> <li>• 3 LED lights</li> <li>• Touch activated and on/off switch</li> <li>• Microfiber lens cloth</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	Always store loupe in included protective case when not in use to avoid damaging the lens.
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Close-up visual inspection of instruments after reprocessing
<b>Preparation for Use</b>	<p><b>Batteries come installed. When changing batteries:</b></p> <ol style="list-style-type: none"> <li>1. Turn item upside down and locate battery compartment. (Please note there are two battery compartments.)</li> <li>2. The compartment closest to the light switch is the one that powers the item.</li> <li>3. The other battery compartment is used to store batteries that are not in use.)</li> <li>4. Pull battery door up and fold door away from the battery compartment <b>Fig. 1</b>.</li> <li>5. Insert two fresh CR2016 coin cell batteries into battery compartment.</li> <li>6. Make sure that the positive side of the battery is facing outward <b>Fig. 2</b>.</li> <li>7. Close battery door securely when finished.</li> </ol>

<b>Diagrams (drawings, pictures):</b>			
	1	2	3
			
	4	5	

<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Select the LED light setting you would like to use (see Fig. 3).             <ol style="list-style-type: none"> <li>a. For constant light, position the switch in the middle position (marked ON).</li> </ol> </li> </ol>
---------------------------------	---

	<p>b. To use the “touch activated” setting, position the switch to the position marked “AUTO”.</p> <p>i. With the switch on AUTO, you simply need to touch the silver ring on top of the loupe that surrounds the lens (as shown in Fig. 4) to activate the light.</p> <p>2. Once you have selected your LED light setting, place the loupe on top of the subject you would like to magnify. Focus the magnifier by rotating the top portion of the loupe until the viewing subject is sharp and clear (see Fig. 5).</p>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• Do not leave batteries installed during extended periods of disuse.</li> <li>• Do not mix new and old batteries. Do not mix alkaline, standard (carbon-zinc) or rechargeable (nickel-cadmium) batteries.</li> <li>• Fire Hazard. Keep magnifier away from direct sunlight to prevent accidental fire.</li> <li>• Do not stare directly at LED lights as this might damage your eyes.</li> </ul>
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	<ul style="list-style-type: none"> <li>• Clean lens with a soft, clean microfiber cloth rubbing gently in a circular motion.</li> <li>• Moisten the lens cloth with water if necessary.</li> <li>• Do not use harsh chemicals, solvents or abrasives as they may damage the lens.</li> </ul>
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	
<b>Reference Documents</b>	ProSys Brochure, ProSys Price List
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com