




**Instructions for Use: Robotic Instrument Guards**

<b>Brand Name of Product</b>	Robotic Instrument Guards
<b>Generic Name of Product</b>	Paper Tip Protectors
<b>Product Code Number(s)</b>	PPT-DV-01, PPT-DV-02
<b>Intended Use</b>	Designed to protect the robotic arm tips during sterilization, transport and storage. Prevents the tips of instruments from penetrating the sterile pouch.
<b>Range of Applications for Product</b>	Robotic instruments
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Single use</li> <li>• Two slot sizes available-             <ul style="list-style-type: none"> <li>○ 15mm</li> <li>○ 30mm</li> </ul> </li> <li>• The wide, smooth design allows for easy insertion and removal from slide a sterilization pouch.</li> <li>• Paper protectors are low linting and have excellent wet and dry mechanical properties to withstand the steam sterilization cycle.</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	100/box, paper
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	To protect the tip of the robotic instrument during sterilization, transport and storage.
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	



<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Fold flap (shorter side) <b>Fig. 1</b></li> <li>2. Toward the base (large size) at the score. <b>Fig. 2</b></li> <li>3. Tuck the flap into the slotted area. <b>Fig. 3</b></li> <li>4. Slide robotic instrument into the first slot, then up through the second slot and back into the pocket formed by the flap and base. <b>Fig. 4</b></li> <li>5. Slide robotic arm with paper tip protector into the paper pouch for sterilization. <b>Fig. 5</b></li> <li>6. Seal pouch as directed by the instrument pouch IFU.</li> </ol>
<b>Interpretation of Test Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	Single-use. Dispose after use in accordance with facility policy.

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• Compatible with steam sterilization.</li> <li>• Recommended cycle is 132°C for 4 minutes, or 135°C for 3 minutes in a pre-vac sterilizer.</li> <li>• Paper Material Compatibility               <ul style="list-style-type: none"> <li>○ Steam up to 137°C</li> </ul> </li> </ul>
<b>Storage</b>	

<b>Additional Information</b>	
<b>Related Healthmark Products</b>	Sterigard Tip Protectors, Vented Instrument Guards, Duo-Vented Guards, Cylindrical Style Instrument Guards
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>