



Instructions for Use: Robotic Instrument Guards

Brand Name of Product	Robotic Instrument Guards
Generic Name of Product	Paper Tip Protectors
Product Code Number(s)	PPT-DV-01, PPT-DV-02
Intended Use	Designed to protect the robotic arm tips during sterilization, transport and storage. Prevents the tips of instruments from penetrating the sterile pouch.
Range of Applications for Product	Robotic instruments
Key Specifications of Product	<ul style="list-style-type: none"> • Single use • Two slot sizes available- <ul style="list-style-type: none"> ○ 15mm ○ 30mm • The wide, smooth design allows for easy insertion and removal from slide a sterilization pouch. • Paper protectors are low linting and have excellent wet and dry mechanical properties to withstand the steam sterilization cycle.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Store in original packaging or in an equivalent clean barrier to maintain cleanliness.
Packaging Contents	100/box, paper
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To protect the tip of the robotic instrument during sterilization, storage and transportation. For use in sterile barrier systems.
Preparation for Use	N/A
Diagrams (drawings, pictures)	

	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Fold flap (shorter side) Fig. 1 2. Toward the base (large size) at the score. Fig. 2 3. Tuck the flap into the slotted area. Fig. 3 4. Slide robotic instrument into the first slot, then up through the second slot and back into the pocket formed by the flap and base. Fig. 4 5. Slide robotic arm with paper tip protector into the paper pouch for sterilization. Fig. 5 6. Seal pouch as directed by the instrument pouch IFU.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	Single-use. Dispose after use in accordance with facility policy.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> • Compatible with steam sterilization. • Recommended cycle is 132 °C for 4 minutes, or 135 °C for 3 minutes in a pre-vac sterilizer. • Paper Material Compatibility • Steam up to 137 °C
Storage	N/A
Additional Information	N/A
Related Healthmark Products	Sterigard Tip Protectors, Vented Instrument Guards, Duo-Vented Guards, Cylindrical Style Instrument Guards

Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	<p>Manufactured By: Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com</p> <p>Distributed By: WiPak Oy Wipaktie 2 P.O. Box 45 15561 Nastola, Finland Tel: +358 20 510 311 Info.fi.wn@wipak.com</p>

2022 Suzanne Latta