


**Instructions for Use: Scissor Test Squares**

<b>Brand Name of Product</b>	Scissor Test Squares
<b>Generic Name of Product</b>	Scissor Sharpness Test
<b>Product Code Number(s)</b>	LT-130, RT608, YT606
<b>Purpose of Product</b>	To test sharpness of scissors
<b>Range of Applications for Product</b>	Scissor sharpness testing
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• LT-130 Latex Free 5 x 5inches– 6.0 gauge</li> <li>• RT608 Latex 216 x 5.5inches– 8.0 gauge</li> <li>• YT606 Latex 216 x 5.5inches – 6.0 gauge</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	<ul style="list-style-type: none"> <li>• LT-130 quantity: 10 Sheets</li> <li>• RT606 quantity: 1 Roll</li> <li>• YT606 quantity: 1 Roll</li> </ul>
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	The test is used to evaluate the sharpness of a scissor.
<b>Preparation for Use</b>	<ul style="list-style-type: none"> <li>• For total length of the scissor that are less than 5 inches, use the yellow test material or the green LT-130 latex-free option.</li> <li>• For total length of the scissor that are greater than 5 inches, use the red material.</li> </ul>
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Make at least two cuts in the test material.</li> <li>2. The scissors should cut all the way through the top of the scissor.</li> <li>3. Start cutting at the middle part of the blade cutting clearly through the tip. If you notice any catching or snagging, tag the scissor for repair/sharpening.</li> <li>4. If the scissors are not sharp enough for use, they should be placed in the appropriate repair bin for future inspection by a qualified instrument repair technician according to your medical facility policy.</li> </ol>
<b>Interpretation of Results</b>	If any instrument does not pass the test, it should be given to a qualified instrument repair technician for a closer inspection.
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• Checking for sharpness in the Central Service Department is part of a quality process.</li> <li>• These are general guidelines to help reduce malfunctioning instruments from being used.</li> <li>• Always remember, if in doubt, place the instrument in the repair bin to have a qualified instrument repair technician check your work.</li> </ul>
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A

<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	Instrument Repair Tags
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>