



**Instructions for Use: Shark Tooth Heat Mitt**

<b>Brand Name of Product</b>	Shark Tooth Heat Mitt
<b>Generic Name of Product</b>	Silicone Heat Mitt
<b>Product Code Number(s)</b>	SLG-001, 09484, 09368
<b>Purpose of Product</b>	For the safe handling of hot trays and racks from autoclaves.
<b>Range of Applications for Product</b>	Trays, racks, and items that are hot being removed from the autoclave.
<b>Key Specifications of Product</b>	<p>Shark Tooth Heat Mitt</p> <ul style="list-style-type: none"> <li>Made of heat deflecting silicone that is sterilizable, washable, and resistant to temperatures from -40 °C (-40 °F) to 230°C (446 °F).</li> <li>The shark tooth ridges provide a safe and secure grip of hot and wet items.</li> <li>Size: SLG-001 Blue: 14 inches length</li> </ul> <p>Silicone Heat Glove:</p> <ul style="list-style-type: none"> <li>Made of silicone, latex free.</li> <li>Protects from surface temperatures ranging from -40- to 464 °F</li> <li>Sizes:             <ul style="list-style-type: none"> <li>09484 Light Blue: 16.90 inches length</li> <li>09368 Light Blue: 13.70 inches length</li> </ul> </li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	Store at room temperature. (Do not expose to direct sunlight)
<b>Packaging Contents</b>	One (1) each per package: <ul style="list-style-type: none"> <li>Shark Tooth Heat Mitt per package</li> <li>Silicone Heat Glove per package</li> </ul>
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	N/A
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	 <p>Shark Tooth Heat Mitt</p>  <p>Silicone Heat Glove</p>
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>Insert your hand completely into the mitt.</li> <li>Mitt or glove can be worn on either right or left hand.</li> <li>Ensure mitt/glove is securely covering your hand for maximum protection.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>Shark Tooth Heat Mitt can withstand temperature of 230 °C (446 °F).</li> </ul>

	<ul style="list-style-type: none"> <li>• Silicone Heat Glove: <ul style="list-style-type: none"> <li>○ The glove can withstand temperature of 240 °C (464 °F).</li> <li>○ Glove will not protect against chemicals.</li> </ul> </li> <li>• Do not clean in an ultrasonic cleaner.</li> </ul>
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	<ul style="list-style-type: none"> <li>• Clean the mitt/glove manually utilizing a neutral pH detergent and a soft bristle nylon brush (e.g., Healthmark brushes 869903N, GDS-1 or GDS-2) and/or a soft non-linting cleaning cloth.</li> <li>• Recommended to use demineralized water for the final rinse.</li> <li>• The mitt may be air dried or dried with a non-linting cloth or towel.</li> </ul>
<b>Cleaning – Automated</b>	<ul style="list-style-type: none"> <li>• Acidic / neutral / alkaline &lt; pH 10.</li> <li>• Consult the detergent manufacturer’s (Mfr.) IFU for the recommended concentration, exposure times, and temperature.</li> <li>• Consult the IFU of the washer Mfr. for correct cycle settings.</li> <li>• Do not clean in an ultrasonic cleaner.</li> </ul>
<b>Disinfection</b>	<ul style="list-style-type: none"> <li>• If desired, after cleaning and drying, the Shark Tooth Heat Mitt may be disinfected utilizing a liquid chemical disinfectant (which is compatible with silicone), in accordance with the chemical Mfr.’s IFUs.</li> <li>• Following disinfection, the mitt should be copiously rinsed with tap water to remove any chemical residuals.</li> </ul>
<b>Drying</b>	The mitt may be air dried or dried with a non-linting cloth.
<b>Maintenance, Inspection, and Testing</b>	<ul style="list-style-type: none"> <li>• Visually inspect mitt gloves for usability and cleanliness.</li> <li>• Damaged mitt/gloves no longer provide protection and must be replaced.</li> </ul>
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<p>If needed, the mitt/glove may be processed through a pre-vac steam sterilizer, in accordance with the healthcare facility’s procedures and conventional cycles at temperatures not to exceed 230 °C (446 °F).</p> <p>Compatible with steam sterilization:</p> <ul style="list-style-type: none"> <li>• 132 °C (270 °F) for 4 minutes</li> <li>• 135 °C (275 °F) for 3 minutes.</li> <li>• Extended cycle: 135 °C (275 °F) for 30 minutes.</li> </ul>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	Insulated Gauntlet, Seamless Autoclave Gloves, Terry Cloth Autoclave Gloves
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>