

Instructions for Use: Shaver Leak Tester

| | |
|--|---|
| Brand Name of Product | Shaver Leak Tester |
| Generic Name of Product | Shaver Leak Tester |
| Product Code Number(s) | SLT-MR, SLT-STK, SLT-STK SMJ, SLT-DYO, SLT-C-L, SLT-ART, SLT-PUMPCS |
| Intended Use | To test the seal of the fluid pathway within arthroscopic shavers. |
| Range of Applications for Product | Checks for performance of interior seal and, if any, leaks are present. |
| Key Specifications of Product | Testing Stop <ul style="list-style-type: none"> • AED-SLT Stryker® Formula - Stop: SLT-STK-SMJ • AED-SLT Midas Rex® - Stop: SLT-MR • AED-SLT Stryker® Formula - Stop: SLT-STK • AED-SLT Dyonics - Stop: SLT-DYO • AED-SLT ConMed® Linatec - Stop: SLT-C-L • AED-SLT-Arthrex® - Stop: SLT-ART • AED-SLT PUMP: SLT-PUMPCS |

| | |
|---|-----|
| Shipping & Storage | |
| Shipping Conditions & Requirements | N/A |
| Storage Conditions | N/A |
| Packaging Contents | N/A |
| Shelf Life | N/A |


| | |
|---------------------------------------|--|
| Instructions for Using Product | |
| Description of Use(s) | To test the seal of the fluid pathway within the arthroscopic shavers. |
| Preparation for Use | N/A |
| Diagrams (drawings, pictures) |  <p>AED-SLT Stryker® Formula SLT-STK-SMJ</p> |



Figure 1




Figure 2




Figure 3



Figure 4




Figure 5




Figure 6

Steps for Use of Product

1. Place the shaver to be tested on a clean flat surface to avoid accidental damage.
2. Select the appropriate testing stop based on the brand of the shaver being tested. Each stop is brand-specific and is not interchangeable (**Fig. 1**)

| | |
|--|--|
| | <p>3. Hold the shaver securely in one hand and the proper testing stop in the other (Fig. 2)</p> <p>4. Insert the testing stop into the distal end of the handpiece until it is fully seated. Align keyed stop(s) with the corresponding slot on the handpiece (Fig. 3)</p> <p>NOTE: Some shaver models may require the engagement of a push-button latch to allow the stop to seat fully.</p> <p>5. Connect the open end of the silicone tubing, attached to the pressure gauge, to the stainless-steel discharge port located at the proximal end of the handpiece. Push together to create an airtight seal (Fig. 4)</p> <p>6. Ensure the pressure relief value of the hand pump is turned fully clockwise to the closed position (Fig. 5)</p> <p>7. Use the bulb of the hand pump to apply pressure:</p> <ol style="list-style-type: none"> No more than 1 to 2 pumps of the bulb will be necessary Observe that the pressure gauge displays and maintains a positive value Positive static pressure on the pressure gauge confirms no leaks in the fluid pathway and considered a passing result (Fig.6) <p>NOTE: Proceed to Step 8 if pressure is not maintained</p> <p>8. Remove testing stop and disconnect silicone tubing from the discharge port of the handpiece</p> <p>NOTE: Inspection, function verification, and sterilization of the shaver can now continue per the manufacturer's IFU.</p> <p>9. If pressure gauge does not maintain static positive value:</p> <ol style="list-style-type: none"> Re-verify all connection points and ensure the pressure relief value is sealed/closed Apply additional pressure using the hand pump Observe the pressure gauge <ol style="list-style-type: none"> If pressure is not holding statically, then shaver is leaking internally The inability to maintain positive pressure is considered a failing result. (Remove all test stops and disconnect tubing from the handpiece) <p>10. Remove failed handpiece from service and identify for repair.</p> |
| Interpretation of Test Results | N/A |
| Contraindications of Test Results | N/A |
| Documentation | N/A |
| Special Warnings and Cautions | <p>Leak testing with Shaver Leak Tester should only be performed after cleaning and before sterilization of the shaver.</p> <p>CAUTION: Please refer to the steps outlined here for proper leak testing techniques. Always review OEM instruction manuals for your device.</p> |
| Disposal | N/A |

| Reprocessing Instructions | |
|--|---|
| Point of Use | N/A |
| Preparation for Decontamination | N/A |
| Disassembly Instructions | N/A |
| Cleaning – Manual | <p>Shaver Leak Testing Stops:</p> <ol style="list-style-type: none"> Rinse excess soil from the stop Using a soft brush, apply detergent (as specified by the manufacturer's IFU) to all surfaces, ensuring complete coverage Rinse thoroughly under running water Apply a small amount of surgical grade lubricant to maintain the O-ring's integrity and help ensure proper seal |

| | |
|---|---|
| Cleaning – Automated | Shaver Leak Testing Stops: <ol style="list-style-type: none"> 1. Load testing stops loosely in a wire basket or similar 2. Select the <i>instrument cycle</i> of the automated washer 3. When unloading, check for the complete removal of any visible soil. If necessary, repeat the cycle or perform manual cleaning 4. A small amount of surgical grade lubricant will maintain the O-ring’s integrity and help ensure proper seal. |
| Disinfection | Shaver Leak Tester Stops and Pump: <ul style="list-style-type: none"> • Remove any visible soil using disposable disinfecting wipe • Thoroughly wipe unit for a minimum of 60 seconds and allow to thoroughly dry before use |
| Drying | N/A |
| Maintenance, Inspection, and Testing | N/A |
| Reassembly Instructions | N/A |
| Packaging | N/A |
| Sterilization | Shaver Leak Testing Stops (as desired) Steam sterilization temperature: <ul style="list-style-type: none"> • 132 °C (270 °F) for 4 minutes, 135 °C (275 °F) for 3 minutes. Maximum temperature 135 °C (275 °F) for 30 minutes NOTE: The hand pump is not compatible with sterilization methods. |
| Storage | N/A |
| Additional Information | N/A |
| Related Healthmark Products | N/A |
| Other Product Support Documents | ProSys™ Brochure, ProSys™ Price List |
| Reference Documents | N/A |
| Customer Service Contact | Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com |