








Instructions for Use: Spray Gun

Brand Name of Product	Spray Gun
Generic Name of Product	Spray Gun
Product Code Number(s)	6020, 80-00-44, 80-10-15, 4240, 4250, 5220, 5110, 5020, 5030, 1HGL4, 4270, 4260, 5070, 4230, 5040, 6000, 6000-6025, 5120, 5010, 6025, 6045, 80-00-013, 6065, 14RVT, 6065H1, 6065H2, 3000, 3005, 5210, 5080, 5085, 5090, 5095, 5100, 5050, 5060, 3070, 4378
Purpose of Product	A tool to use in the cleaning of items with water or air.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> • Use with the “universal” slip-tip, or the optional quick connect set of tips for specialized functions. • Maximum water or air pressure six (6) bar (86 psi).
	N/A

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Conditions	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	A tool to use in the cleaning of items.
Preparation for Use	<ul style="list-style-type: none"> • Best installed by a licensed plumber. Please contact Facilities and/or Engineering for installation of the Spray Gun. Installation of the Spray Gun is performed by the healthcare facility. • Connection requirements <ul style="list-style-type: none"> ○ Water connection: <ul style="list-style-type: none"> ▪ Plumbing threaded tap or valve: 3/8-, 1/2-, or 3/4-inch ▪ Standard gun kit includes 3/4-inch female connection. ○ Air connection: <ul style="list-style-type: none"> ▪ If operated with compressed air, the tubing must be connected to the air supply by a suitable connector (Note: cannot be supplied by Healthmark). ▪ When connecting the gun to water or air below the counter, it is recommended that a rosette be installed. It will minimize potential friction damage to the hose. A 1 3/8 inches diameter hole will need to be drilled in the counter or sink top to accommodate the rosette. • Prior to operation, the user must make sure attachments are correctly seated on the safety cone. • If the Spray Gun is not in use, ensure that water or air pressure have been turned off. • When the hose of the Spray Gun is installed to the water line, a separate shutoff valve needs to be installed to the spigot. (The Spray Gun is not a controller.)
Diagrams (drawings, pictures)	

		
	<p><i>Figure 1: 6000-6025 Quick Connect Tips</i></p>	<p><i>Figure 2: 5220 Specialized Tips</i></p>
		
	<p><i>Figure 3: 3000, Nylon Spray Gun Wall Mount</i> <i>Figure 4: 3005, Chrome Brass Spray Gun Wall Mount</i></p>	<p><i>Figure 5: 5210 Spray Gun Brush Set- 7 Brush Attachments and Holder</i></p>
		
	<p><i>Figure 6: Threaded nut for controlling maximum flow</i></p>	<p><i>Figure 7: Flow rate controlled by trigger handle</i></p>
		<p><u>Key to Spray Gun Components</u></p> <ol style="list-style-type: none"> 1. Gun Rinser Tip Holder complete kit 2. "O" Ring black 3 mm x 1 mm 3. Spiral Tension Pin 4. Trigger Handle 5. Grip for Spray Gun 6. Hose Connector 7. Gun Cap 8. Flow Control Kit
	<p><i>Figure 8: Diagram Spray Gun Construction</i></p>	
<p>Steps for Use of Product</p>	<p>Prior to operation, the user must:</p> <ol style="list-style-type: none"> 1. Take care not to point Spray Gun at self/others orifices prior to operation. 2. Make sure attachments are correctly seated on the safety cone. <ol style="list-style-type: none"> a. Firmly push the attachments onto the Spray Gun's rinse tip. b. Ensure the valve is open that is attached to the hose. <p>When in operation:</p> <ol style="list-style-type: none"> 3. Wear appropriate PPE per the facility's policy. 4. Set Spray Gun's maximum pressure with small threaded nut located before the trigger handle. 	
<p>Interpretation of Results</p>	<p>N/A</p>	
<p>Contraindications of Test Results</p>	<p>N/A</p>	
<p>Documentation</p>	<p>N/A</p>	
<p>Special Warnings and Cautions</p>	<p>Cautions:</p>	

	<ul style="list-style-type: none"> • Verify the safe and effective performance of any tubing chosen for use other than that supplied by Healthmark. The liability for damage or injury is the responsibility of the facility. • Warranty claims will <i>only</i> be honored for components purchased from Healthmark. • Maximum water or air pressure to use is 6 bar (86 psi). • Be sure to wear proper PPE, according to facility policy. <p>Warnings:</p> <ul style="list-style-type: none"> • Do not use Spray Gun for any medical treatments. • Do not point the Spray Gun at parts or orifices of the body when in use. Compressed water or air jet may be harmful. •
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	
Disassembly Instructions	N/A
Cleaning – Manual	<ul style="list-style-type: none"> • Use only warm water and a mild detergent to clean the Spray Gun with all add-ons. • The user should not operate the Spray Gun with demineralized water since it causes corrosion.
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	<ul style="list-style-type: none"> • It is recommended to decalcify the Spray Gun approximately every four (4) weeks. This is beneficial for the service life and reliability of the Spray Gun. <ul style="list-style-type: none"> ○ The gun should be disassembled. ○ Inner parts are placed in a decalcifying agent. • It is not necessary to carry out water quality tests in the case of using potable water. • Replacement parts and instructions for repair are available from Healthmark.
Reassembly Instructions	<ul style="list-style-type: none"> • The complete reassembly of the clean gun must be carried out under sterile conditions. • Before screwing the individual parts together, make sure that all O-Rings are in perfect condition. • Sealing elements must be moistened with a suitable physiological saline solution before use. • Sealing rings are to be mounted at the intended locations and locked in place by pressing them firmly by hand.
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com