



	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Remove protective strip.</li> <li>2. If desired, Place indicator inside pouch (Type1 Indicator).</li> <li>3. Fold flap.</li> <li>4. Press down from the center with both thumbs and move outward. <ol style="list-style-type: none"> <li>a. Apply even pressure along the entire length of the strip (including the entire length of the outer folded edge), while moving thumbs.</li> <li>b. Repeat at least three times more to ensure a proper seal.</li> </ol> </li> <li>5. Examine the seal to make sure you do not see any uneven spots or ridges.</li> <li>6. Check each end of the seal to be sure no channel exists, and the seal is complete.</li> <li>7. Apply pressure again (evenly) to smooth out uneven spots, ridges, or channels.</li> <li>8. Peel pouches should be located in the sterilizer on a rack and stood vertically on their edges to maximize exposure to the sterilization process.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	<ul style="list-style-type: none"> <li>• The Steriking® see-through range of peel packages conform to the international product standards and norms: ISO 11607-1, ISO 11607-2.</li> <li>• The products are registered by FDA under 510(k) Premarket Submission Nos.: K803293, K953776 and K973827. Wipak Oy™ is certified to ISO 9001:2000; ISO 14001: 2004; OHSAS 18001: 1999; ISO 22000: 2005 and DS 3027: 2002.</li> <li>• Steriking® sterilization packages are designed, validated, and manufactured to suit their intended purposes.</li> </ul>
<b>Special Warnings and Cautions</b>	The Steriking® standard range of see-through packaging is not suitable for sterilization by irradiation or by hot, dry air at the temperatures over 140 °C (284 °F).
<b>Disposal</b>	Please refer to the local/national regulations regarding waste disposal.

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instruction</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	Be sure to check the seal is complete from end-to-end.
<b>Reassembly Instruction</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• Compatible with steam and EO sterilization.</li> <li>• <i>Pre-vacuum Sterilizer:</i> Wrapped trays (and instruments) should be exposed to: <ul style="list-style-type: none"> <li>○ 132 °C (270 °F) for 4 minutes</li> <li>○ 135 °C (275 °F) for 3 minutes.</li> </ul> </li> <li>• Compatible with Gravity Steam Sterilization at 121 °C (250 °F) for 30 minutes.</li> <li>• Arrange pouches so that there is minimal to no contact between pouches.</li> <li>• Chemical Indicators conform to ISO 11140-1 Class 1: Process indicators. <ul style="list-style-type: none"> <li>○ Steam indicator changes color from blue to dark brown/black (New Design from pink to brown)</li> <li>○ EO gas indicator from pink to yellow/orange on the outside of the pouch (New Design from purple-pink to beige-orange). Color can vary depending upon the concentration of EO gas and amount of humidity.</li> </ul> </li> </ul>
<b>Storage for Sterilized Pouches</b>	<ul style="list-style-type: none"> <li>• Recommended long-term storage conditions: <ul style="list-style-type: none"> <li>○ Temperature: 18- to 24 °C (64.40- to 73.40 °F)</li> <li>○ Humidity 30–60% rh</li> <li>○ Do not exceed 40 °C for sterile storage area.</li> </ul> </li> </ul>

	<ul style="list-style-type: none"> <li>Facilities should develop a policy identifying events that may compromise the sterility of packaged items. This could include a maximum period of time for storage.</li> <li>Sterilized self-seal pouches were tested to 3 years of accelerated and real time aging.</li> </ul>
<b>Additional Information</b>	<ul style="list-style-type: none"> <li>These products are for single use only.</li> <li>Maximum tested sterilization time: 30 minutes.</li> <li>Maximum tested temperature: 135 °C.</li> <li>Validation testing are performed for single and double pouching.</li> </ul>

<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	Marking and labeling Steriking® Packaging- <a href="https://www.hmark.com/wp-content/uploads/2022/03/Marking-and-labeling-Sterking-Packaging.pdf">https://www.hmark.com/wp-content/uploads/2022/03/Marking-and-labeling-Sterking-Packaging.pdf</a>
<b>Reference Documents</b>	Sterilization Packaging Brochure, Sterilization Packaging Price List
<b>Customer Service Contact</b>	US Distributor: Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com