

For more information about ProFormance™ products and verification of the cleaning process, call 800-521-6224 or visit HMARK.com

Make Clean Seen™

Instruction Guide for Ultrasonic Cleaners, Scope Washers



How to Use this Booklet

- 1 Determine the test or test kit you are using.
- 2 If you are using a kit, identify each individual test in the kit.
- 3 Read and follow the instructions for each test to be used.
- 4 Perform equipment inspection as specified on the log sheet.
- 5 Record results on the log sheet.
- 6 Immediately report any test failure to department management.


1

AquaTest™

Reorder# AT101

Follow these directions first for Cold, then Hot Water:

- 1 Dip entire strip into water for 5 seconds, then remove.
- 2 Shake once briskly to remove excess water from the test strip.
- 3 Wait 20 seconds.
- 4 Compare color within 10 seconds to pH, TOTAL ALKALINITY, TOTAL HARDNESS on the interpretation chart.
- 5 Record results on log sheet.
- 6 Report any deviations from expected values.



2

Temp-LC2

Reorder#
TEMP-LC2

- 1 Use one (1) Temp-LC2 per ultrasonic.
- 2 Peel thermometer from release paper.
- 3 Apply to the inside, side wall of the tank below the intended fluid line and where the thermometer can be easily read.
- 4 Temp-LC2 should be used to adjust the water temperature to the proper, optimal range for the cleaning agent used as specified by the manufacturer.
- 5 A color of green indicates the temperature is at the exact value. Blue is the temperature just below and tan is the temperature just above the actual temp.
- 6 Monitor the temperature during use of sonic to insure compliance with target temperatures. Immediately report failure to achieve the target.
- 7 Record observations on log sheet.



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

PB-04/22

healthmark www.hmark.com
°F 60 70 80
°C 16 21

3

SonoCheck™

Reorder# TI108

- 1 Be sure to run a degassing cycle prior to testing.
- 2 Depending upon the size of the sonic bath, arrange up to 3 SonoChecks in the appropriate loading basket (See chart below).



- 3 Program the sonic for a standard instrument cleaning cycle (follow manufacturer instructions).
- 4 At end of cycle, check the SonoChecks for the appropriate color change.
- 5 Record results on log sheet.
- 6 Report if color indicates a failed test.


| | | |
|---|---|---|
| | | X |
| | X | X |
| X | X | X |

4

TOSI®

Reorder# WT101


- 1 Be sure to run a degassing cycle prior to testing.
- 2 Place one (1) TOSI® test in the center of the appropriate loading basket and load according to instructions.
- 3 Program the sonic for a standard instrument cleaning cycle (follow manufacturer instructions).
- 4 At the end of the cycle, visually examine the test for any remaining residual soil and compare to the Interpretation Guide.
- 5 Record results on log sheet.
- 6 If any residual soil remains, (score of 1 - 5 on the Interpretation Guide) immediately report results.



5

LumCheck™

Reorder# WLC-101




- 1 Unscrew and open the LumCheck™ device. Detach part A from part B.
- 2 Open the foil pouch and insert the test object part C into part B.
- 3 Screw part A to part B close.
- 4 Connect the LumCheck™ device to one of the washer's channel irrigators and process according to the manufacturer's instructions.
- 5 Load the lumen washer/disinfector normally.
- 6 Process as per normal procedure.
- 7 After the process is complete, disconnect the LumCheck™ and examine the test for any remaining residual soil and compare to the Interpretation Guide. Report any failure.
- 8 Record results on log sheet.



6

FlexiCheck™

Reorder# FLEXCHK-101



- 1 Unscrew FlexiCheck™ – device; detach part A from part B.
- 2 Open foil pouch of FlexiCheck™ and insert the test object part C into part B. Do not touch the test soil.
- 3 Screw part A to part B.
- 4 Connect the FlexiCheck™ device to one of the washer's channel irrigators and process according to the manufacturer's instructions.
- 5 After the process is complete, disconnect part A from part B and remove part C, the coupon without touching the soil area.
- 6 Examine the test for any remaining residual soil and compare to the Interpretation Guide. Report any failure.
- 7 Record results on log sheet.

