

Brand Name of Product	Gorilla Bags™
Generic Name of Product	Dust Covers
Product Code Number(s)	GB-001 PP, GB-002 PP, GB-003 PP
Intended Use	To be used as protective packaging (as defined by ISO 11607-1)
Range of Applications for Product	Wrapped trays
Key Specifications of Product	Manufactured from spun polypropylene.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Clean, dust-free area, and away from fluorescent or ultraviolet light.
Packaging Contents	N/A
Shelf Life	Five (5)-year shelf life prior to use when stored at: <ul style="list-style-type: none"> • Humidity: 20–80% rh • Temperature: 10 °C (50 °F) to 35 °C (95 °F).

Instructions for Using Product	
Description of Use(s)	Added protection from environmental conditions as protective packaging.
Preparation for Use	Select proper size bag for wrapped tray.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Wrap tray (according to policy) with facility approved wrap. 2. Once the tray is wrapped, place bag on table. 3. Place wrapped tray into Gorilla Bag™. 4. Fold the opening of the bag forward three times. 5. Seal with two strips of autoclave tape <ol style="list-style-type: none"> a. Place each strip on the bag b. Wrap tape underneath the fold to secure the closure of the bag. 6. After cool down, inspect the Gorilla Bag™ before it is removed from the rack. 7. Follow facility polic(ies) with proper handling techniques by lifting the wrapped tray off the rack and set onto the transportation cart to be stored. (Note: Do not drag the tray.) 8. After the set has been transported from storage to the OR, remove Gorilla Bag™.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	<ul style="list-style-type: none"> • Is not a microbiological layer. • Never grab/hold the end of the Gorilla Bags™ while transporting/carrying items. This action can compromise the integrity of the bag. • Proper care and handling of the Gorilla bags™ should be observed at all times.
Disposal	<ul style="list-style-type: none"> • Single-use only. • Landfill or incinerate based upon state or local regulations. • Gorilla Bags™ are composed of polypropylene plastic, which has a recycling code 5 ♻️.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A

Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Gorilla Bags™ are completely porous to the sterilant and can be used in pre-vacuum steam sterilizers for a maximum temperature of 138 °C (280.4 °F) for a maximum time of 30 minutes.
Storage	N/A
Additional Information	Available in four (4) different colors: purple, red, orange, and yellow.

Related Healthmark Products	Sterilization Packaging Brochure, Sterilization Packaging Price List
Other Product Support Documents	N/A
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com