



**Instructions for Use: Paper Tip Protectors**

<b>Brand Name of Product</b>	Paper Tip Protectors
<b>Generic Name of Product</b>	Paper Point, Plastic & Needle Protectors
<b>Product Code Number(s)</b>	PP-001, PP-002, PP-003, PP-004, PP-005, PP-006, PO-007C, PO-008C
<b>Purpose of Product</b>	Designed to protect instrument tips during sterilization, transport and storage.
<b>Range of Applications for Product</b>	Slot Diameter ranges from 5mm to 30mm to accommodate a variety of sizes.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Single use.</li> <li>• Wide, smooth design allows for ease to slide in and out of instrument pouches.</li> <li>• Paper protectors are low linting and have excellent wet and dry mechanical properties to withstand the steam sterilization cycle.</li> <li>• Plastic protectors are see-through for the protection of tips and points where visual identification is needed.</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	<ul style="list-style-type: none"> <li>• 200/box, plastic: (PO-007C, PO-008C)</li> <li>• 1000/box, paper: PP-002, PP-004, PP-006</li> <li>• 2000/box, paper: PP-001, PP-003, PP-005</li> </ul>
<b>Shelf Life</b>	Five (5)-year shelf life.

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	To protect sharp instruments during sterilization, transport, and storage.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	N/A
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Fold flap (shorter side) toward the base (large size) at the score.</li> <li>2. Tuck the flap into the slotted area.</li> <li>3. Slide instrument into the first slot, then up through the second slot, and finally back into the pocket formed by the flap and base.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	<ul style="list-style-type: none"> <li>• Single use.</li> <li>• Dispose after use in accordance with laws and facility regulations.</li> </ul>

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• Compatible with steam sterilization.</li> <li>• Recommended cycle in a pre-vac sterilizer: <ul style="list-style-type: none"> <li>○ 132 °C (270 °F) for four (4) minutes</li> <li>○ 135 °C (275 °F) for three (3) minutes.</li> </ul> </li> <li>• Paper Material Compatibility: <ul style="list-style-type: none"> <li>○ Steam up to 137 °C (279 °F)</li> </ul> </li> </ul>

	<ul style="list-style-type: none"> <li>○ Ethylene Oxide.</li> <li>● Plastic Material Compatibility: <ul style="list-style-type: none"> <li>○ Steam up to 137 °C (279 °F)</li> <li>○ Ethylene Oxide</li> <li>○ VH2O2</li> </ul> </li> </ul>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	Translucent Tip Protectors, Vented Instrument Guards, Duo-Vented Guards, Cylindrical Style Instrument Guards
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>