








**Instructions for Use: Weinstein Instrument Holder**

<b>Brand Name of Product</b>	Weinstein Instrument Holder
<b>Generic Name of Product</b>	Instrument Stringers
<b>Product Code Number(s)</b>	HE1-0010-W, HE1-0012-W, HE1-0014-W, TBIS-8
<b>Purpose of Product</b>	This product is used to expand instruments to better expose hard-to-reach areas during cleaning and sterilization.
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Manufactured of stainless steel.</li> <li>Available in a variety of widths and lengths.</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	Makes hard to reach areas easier to clean and sterilize. Keeps instruments maximally exposed for hard-to-reach areas during cleaning and sterilization.
<b>Preparation for Use</b>	N/A

<b>Diagrams (drawings, pictures)</b>	
	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  <p>1</p> </div> <div style="text-align: center;">  <p>2</p> </div> <div style="text-align: center;">  <p>3</p> </div> </div>
	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  <p>4</p> </div> <div style="text-align: center;">  <p>5</p> </div> </div>

<b>Steps for Use of Product</b>	1. Place one hand at each end of the stringer.
---------------------------------	--

	<ol style="list-style-type: none"> <li>2. Gently squeeze the stringer arms inwards and pull to separate the stringer from the top.</li> <li>3. String handled instruments (in set assembly order) onto the stringer by placing one ring on each of the stringer arms. <ol style="list-style-type: none"> <li>a. Continue until all instruments are loaded onto the stringer.</li> <li>b. Lock the stringer by reattaching the top. (NOTE: You may need to gently squeeze one of the arms inward to thread the crescent end.)</li> </ol> </li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	<ul style="list-style-type: none"> <li>● Before initial use of the stringer, cleaning is required before initial use. <ol style="list-style-type: none"> <li>1. Manually wash using neutral pH detergent and warm water.</li> <li>2. Place the stringer(s) in ultrasonic cleaner for effective cleaning.</li> <li>3. Next, place stringer(s) in automated washer system for effective cleaning.</li> </ol> </li> <li>● After initial use, the stringer is adequately cleaned during the cleaning of the instruments.</li> </ul>
<b>Cleaning – Automated</b>	See instructions under “Steps for Use of Product” for directions on how to string instruments.
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ol style="list-style-type: none"> <li>1. See instructions under “Steps for Use of Product” for directions on how to string instruments.</li> <li>2. When using for sterilization of instruments: <ol style="list-style-type: none"> <li>a. Open to the optimal width for use in the instrument set.</li> <li>b. Consult the instrument manufacturer’s IFU for specific instructions.</li> </ol> </li> <li>3. Stringers are materially compatible with steam and EO sterilization.</li> <li>4. Pre-Vac sterilization: <ol style="list-style-type: none"> <li>a. Four (4) minutes at 132 °C (270 °F)</li> <li>b. Three (3) minutes at 135 °C (275 °F).</li> </ol> </li> <li>5. Gravity steam sterilization: 30 minutes at 121 °C (250 °F).</li> <li>6. Maximum steam sterilization at 135 °C (275 °F) for 30 minutes.</li> </ol>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	Stringers
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com