





Instructions for Use: FIS Wall Mount with Bracket

Brand Name of Product	FIS Wall Mount with Bracket
Generic Name of Product	Wall Mount with Bracket
Product Code Number(s)	FIS-WMB
Intended Use	To hold the FIS Tablet Case while securely mounted to the wall.
Range of Applications for Product	To hold FIS Tablet Case
Key Specifications of Product	<ul style="list-style-type: none"> • Made of stainless steel with powder coat finish. • +/- 10° for optimal viewing angle. • Compatible mounting points- 75 mm & 100 mm • Length- 4.75 inches (121 mm) • Width- 1.65 inches (42 mm) • Weight- 14/6 ounces (414 g) • Maximum load capacity- 30 pounds (13/6 kg).

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	<ul style="list-style-type: none"> • One (1)-Plate • One (1)-Tilting Wall Plate • Four (4) Drywall Anchors • M5 Inner Hexagon Spanner • Four (4) Phillips Pan Head Screws. • Four (4) M4 x 0.7- x 12 mm • Four (4) M4 x 0.7- x 20 mm Short Pan Head Screws • Four (4) Washers 6.2- x 16- x 1.5 mm
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Holds the FIS Tablet Case while securely mounted to the wall.
Preparation for Use	<ol style="list-style-type: none"> 1. Remove parts from packaging. 2. A Phillips screwdriver is required as assemble the FIS Wall Mount Bracket.
Diagrams (drawings, pictures)	<div style="display: flex; flex-wrap: wrap;"> <div style="width: 50%;">  <p style="text-align: center;">Figure 1</p> </div> <div style="width: 50%;">  <p style="text-align: center;">Figure 2</p> </div> <div style="width: 50%;">  <p style="text-align: center;">Figure 3</p> </div> <div style="width: 50%;">  <p style="text-align: center;">Figure 4</p> </div> </div>

	<p style="text-align: center;">Figure 5</p>
Steps for Use of Product	<ol style="list-style-type: none"> 1. Drill holes where the wall mount will be placed with the 4 drywall anchors. (Fig.1). (Note: Use a level to ensure the tilting wall plate is horizontal. 2. Ensure the plate is securely mounted to the wall. (Fig. 2). 3. With the plate securely mounted, hang the entire assembly onto the tilting plate. (Fig. 3). 4. Make sure the 2 locking tabs on the bottom plus the 1-screw hole on top of the plate all lined up. (Fig. 3). 5. After the plate is installed and hanging on the wall plate, use the M4 x 12 mm Phillips pan head screws to secure these two parts. (Fig. 4). 6. Adjust the tilt angle so that the tablet enclosure is viewable in its ideal state. Maximum angle is -5° to $+15^{\circ}$. (Note: The plate can be adjusted for optimal viewing angle. Once set to the desired angle, the mount can be easily locked by tightening one of the cap screws.) (Fig. 5). 7. After adjusting the tilt/tablet viewing angle, use the included M5 inner hexagon spanner to secure the hexagon screw. (Fig. 5).
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Wipe down with isopropyl alcohol.
Cleaning – Automated	N/A
Disinfection	Use product compatible with stainless steel.
Drying	N/A
Maintenance, Inspection, and Testing	Tighten or adjust, as necessary, for optimal operation and safety.
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A
Related Healthmark Products	FIS Mobile Cart, Flexible Inspection Scope FIS-007
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com