


## 1 AquaTest™

Reorder# AT101


Follow these directions first for Cold Water, then Hot Water:

1. Dip entire strip into water for 5 seconds, then remove.
2. Shake once briskly to remove excess water from test strip.
3. Wait 20 seconds.
4. Compare color within 10 seconds to pH, Total Alkalinity, Total Hardness on the interpretation chart.
5. Record results on the log sheet.
6. Report any deviations from expected values.



## 4 LumCheck™/FlexiCheck™


1. Unscrew and open the LumCheck™ or FlexiCheck™ device. Detach part A from part B.
2. Insert the test object part C into part B.
3. Screw part A to part B of the LumCheck™ or FlexiCheck™ device closed.
4. Connect the LumCheck™ or FlexiCheck™ device to the rinse system of the washer and process the same as instruments.
5. Load the washer/disinfector normally.
6. Process as per normal procedure.
7. After the cleaning process is completed, disconnect the LumCheck™ or FlexiCheck™ and examine the test object for visual cleanliness.
8. Record results on log sheet.



Reorder: LumCheck# WLC-101 / FlexiCheck# FLEXCHK-101

## 2 TempaChek™-90

Reorder# TEMP-90




1. Use one (1) TempaChek™-90 per washer.
2. Peel thermometer from release paper.
3. Apply to any clean surface ensuring that the indicator has adhered to the surface.
4. TempaChek™-90 should be read immediately after the Cold Water Rinse stage.
5. Record results on log sheet.
6. Temperature should not exceed 100°F. Immediately report any result that exceeds this temperature.

## 5 SonoCheck™

Reorder# TI108




- If the washer has an ultrasonic bath...
1. Confirm that the ultrasonic runs a degas cycle.
  2. Depending on size of sonic bath, arrange 3 SonoChecks in the appropriate loading basket (see chart to the right).
  3. Run the washer with ultrasonic through an instrument cycle.
  4. At end of the cycle, compare SonoCheck for the appropriate color change.
  5. Record results on log sheet.
  6. Report if color indicates a failed test.



## 3 TempaChek™-170

Reorder# TEMP-170



1. Use one (1) TempaChek™-170 on each level of the instrument rack.
2. Peel thermometer from release paper.
3. Apply to any clean surface ensuring that the indicator has adhered to the surface.
4. TempaChek™-170 should be read after the Thermal Disinfection stage and before the Drying stage.
5. Record results in log sheet.
6. Report any deviation from targeted temperature on interpretation guide.

## 6 TOSI®

Reorder# WT101

1. Secure one (1) TOSI® per level of the instrument rack in the center of an empty tray (position "A" in the diagram to right).
2. Process using the "instrument" wash cycle.
3. If using with TempaChek™-170, you may remove test prior to the Drying Cycle.
4. Compare the test to the Interpretation Guide.
5. Record results on log sheet.
6. If any residual soil remains, (score of 1-5 on the Interpretation Guide), immediately report results.

