




<b>Brand Name of Product</b>	Magnetic Field Tester
<b>Generic Name of Product</b>	Magnetic Field Tester
<b>Product Code Number(s)</b>	MFT-KIT, MFT-001, MFT-CASE, MFT-FILM
<b>Intended Use</b>	Used to detect the presence of magnetization on medical instruments.
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	<b>MFT-KIT</b> – Magnetic Field Tester kit includes: 2 MFT-CASE and 6 MFT-FILM <b>MFT-001</b> – Magnetic Field Tester includes: 1-MFT-CASE and 1-MFT-FILM <b>MFT-CASE</b> – Case that snaps together around MFT-FILM <b>MFT-FILM</b> – Magnetic film strip.

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	Keep away from sunlight, heat, and radioactive sources.
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	To detect the presence of magnetization on medical instruments.
<b>Preparation for Use</b>	<p>The housing consists of two (2) plastic pieces that snap together. Between the two (2) pieces is a magnetic film strip that will detect magnetization.</p> <ol style="list-style-type: none"> <li>1. Place the instrument to be inspected on the black side of the film while you are looking at the green side of the film.</li> <li>2. If the instrument is magnetized, you will see a color change (like a rainbow) on the green side.</li> <li>3. It is recommended to demagnetize the instrument prior to sterilization.</li> </ol>
<b>Diagrams (drawings, pictures)</b>	<div style="display: flex; flex-wrap: wrap; justify-content: space-around;"> <div style="text-align: center;">  <p>Figure 1</p> </div> <div style="text-align: center;">  <p>Figure 2</p> </div> <div style="text-align: center;">  <p>Figure 3</p> </div> </div> <p style="text-align: center; margin-top: 10px;">Magnetic Field Tester</p>




Figure 4




Figure 5



Figure 6




Figure 7




Figure 8




Figure 9




Figure 10




Figure 11




Figure 12

#### Steps for Use of Product

1. Hold the Magnetic Field Tester in one hand and the instrument you would like to check in the other. (**Fig. 1**).
2. *Gently* glide the instrument over the back side (black side) of the tester. (**Figs. 2–8**).
3. Observe the green side of the film for test results.
  - Fail = color change of the green side of the film. (**Fig. 9 and 10**).
  - Pass = no color change on the green side of the film. (**Fig.11, 12**).
4. If the instrument fails the magnetization test, the instrument should be demagnetized prior to sterilization using the department's demagnetizer.
5. After demagnetization, check the instrument again with the Magnetic Field Tester to make sure it has been properly demagnetized. (Pass).

#### Interpretation of Test Results

N/A

#### Contraindications of Test Results

N/A

<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• DO not sterilize the Magnetic Field Tester.</li> <li>• Do not use if the package is damaged.</li> <li>• The magnetic film strip should be replaced when damage becomes visible and obvious, or every three (3) months.</li> <li>• The housing should be used until it is no longer functional as a film holder.</li> </ul>
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	The housing and film can be cleaned using an IPA Wipe at the end of each shift or as needed.
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	ProSys™ Brochure
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark, A Getinge company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com