




Instructions for Use: Weisenbach Instrument Forceps Stringer

Brand Name of Product	Weisenbach Instrument Forceps Stringer
Generic Name of Product	Stringers
Product Code Number(s)	10-1775, 10-1776
Intended Use	Keeps ring handled instruments maximally exposed during sterilization.
Range of Applications for Product	Used for ring handled instruments for sterilization.
Key Specifications of Product	Manufactured of stainless steel.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Keeps ring handled instruments maximally exposed during sterilization.
Preparation for Use	<ol style="list-style-type: none"> 1. This requires two stringers to hold the rings open during sterilization. 2. Each bar will hold one ring handle for that instrument.
Diagrams (drawings, pictures)	 <p style="text-align: center;">Figure 1</p>  <p style="text-align: center;">Figure 2</p>

	 <p style="text-align: center;">Figure 3</p>
Steps for Use of Product	<ol style="list-style-type: none"> 1. Unlock the “T” of the stringer by placing it in a vertical position. Fig. 1 2. Place the ring handled instruments in the open position. 3. The rings of the instrument are put through both bars of the stringer. Fig. 2 4. When finished placing the instruments through the bars, the “T” is then locked back in the horizontal position. Fig. 3
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Before initial use of the stringer, cleaning is required. Manually wash stringer(s) with a neutral pH detergent and warm water.
Cleaning – Automated	Use the “instrument cycle” of the washer-disinfector compliant to ISO 15883-1 & 2.
Disinfection	The thermal disinfection stage of a washer-disinfector compliant to ISO 15883-1 & 2.
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> • See instructions under “Steps for Use of Product” for directions on how to string instruments. • When using for sterilization of instruments, open to the optimal width for use in the instrument set. Consult the instrument manufacturer’s IFU for specific instructions. • Stringers are materially compatible with steam and EO sterilization. • Steam sterilization at 132°C (270°F) for 4 minutes or 135°C (275°F) for 3 minutes. • Maximum steam sterilization at 135°C (275°F) for 30 minutes.
Storage	N/A
Additional Information	N/A
Related Healthmark Products	Stringers
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

