

Triple Enzyme Detergent

Sanizyme

Sanizyme is our low-foaming, proprietary blend of enzymatic detergents that removes proteins and debris from surgical instruments. This special blend of three enzymes is pH-neutral, making it highly effective in removing bioburden.



Indications

Sanizyme is a low-foaming, triple enzyme detergent that is pH-neutral, making it highly effective in removing organic material safely and effectively from surgical instruments. This product is suitable for ultrasonic cleaners.

Instructions for Use

1. Fill the ultrasonic bath to an appropriate level with warm water. For optimal results, wash temperature should be set at 95-105°F (35-40°C), or in accordance with your equipment manufacturer IFU, ensuring the quantity of water used allows the correct dose to be calculated. Dosing is 1 oz./gal (8 ml/l), depending on soil load and water quality. Higher doses will be required for higher soil loads and for poor water quality.
2. Run the ultrasonic cleaning cycle, ensuring sufficient time to obtain the desired cleaning. Longer cleaning times may be required for heavily soiled items.
3. Discard the solution in accordance with facility policy. Solutions should not be reused if they are visibly soiled.

Recommended Use

For removing organic material found on surgical instruments, such as blood, tissue, mucus, feces, and other organic material.

Storage and Handling

- Sanizyme must be stored in a cool, dry place. Keep partially used containers closed and away from heat and humidity.
- Optimum storage temperatures 33°-77°F (1°-25°C).
- Store upright and in original container.
- Use with adequate ventilation. Avoid breathing vapors, spray, or mist.
- Keep out of reach of children

Size	Part number
2.5 gallons (9.46 liters)	6003000833
1 gallon (4-pack, 15.16 liters)	6003000834
5 gallons (18.91 liters)	6003000835
15 gallons (56.78 liters)	6003000836

Cautions

Warning: May cause skin and eye irritation. Wash hands thoroughly after handling. Wear eye protection, protective gloves. IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously with water for several minutes.

- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.
- Consult the manufacturer's recommended instructions to determine additional processing steps necessary.
- Dispose of contents and container in accordance with local, regional, national, and international regulations.

First Aid:

Eye Contact: In case of contact, immediately flush eyes with plenty of water. Remove contact lenses and flush again. Get medical attention if irritation persists.

Skin Contact: In case of contact, immediately flush skin with soap and water. Wash clothing before reuse.

Inhalation: If inhaled, remove to fresh air. If exposed person is not breathing, give artificial respiration or oxygen applied by trained personnel. Get medical attention.

Ingestion: Do not induce vomiting. Never give anything by mouth to an unconscious person. If irritation persists, get medical attention. Dispose of contents and container in accordance with local, regional, national, and international regulations.

Do not mix with other chemicals.

Store in a cool, dry place.

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries.

Copyright 2026 Getinge AB or its subsidiaries or affiliates. All rights reserved.

Manufacturer: Vintage Chemical, Inc, 2007 Bremer Rd, Fort Wayne, IN 46803
Distributed By: Healthmark Industries, Co, 18600 Malyn Blvd, Fraser, MI 48026

FILE-004963-V1