



## **FOR IMMEDIATE RELEASE:**

### **CONTACTS:**

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **12 Day Indicator**

Fraser, MI March 7, 2019 - Healthmark Industries is excited to announce the 12 Day Indicator Hangtime Label to its Endoscopy product line.

Designed as a visual reminder for healthcare workers, the label features a twelve-day elapsed time indicator adhered to a blue 8.5 x 2.5-inch self-looping hangtime label, allowing for quick and ongoing observable clarification of how long it's been since an endoscope has been reprocessed. The label allows staff to document critical information including the scope serial number and healthcare worker initials, as well as the date and time the indicator was activated.

No preconditioning is required, there are no moving parts, and there is no power source. To activate the indicator, firmly press the blister on the strip, and upon activation a safe dye migrates along a twelve-day display run-out window, which shows elapsed time (Accurate to within +/- 10%).

Additionally, Healthmark Industries offers the related products, including the 1 Hour and 7 Day Indicator Hangtime Labels: "HTK-1H" (*One Hour Indicator*) and "HTK-7D2" (*7 Day Indicator Hangtime Label*).

Visit [www.hmark.com](http://www.hmark.com) or call 800-521-6224 for more information.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.