



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Offers a New Product to Combine with Artificial Test Soil

Fraser, MI January 28, 2014 - Healthmark Industries has announced that they are now offering a product called, "BactiQuik." It is a new prepackaged 1 ML vial of microbial suspensions that can be quickly and easily mixed with "Artificial Test Soil" (ATS) to perform validation testing for disinfection and sterilization.

The vials that they come in contain a pure population of approximately 1×10^8 CFU/ML. They are offered individually by either viable *Geobacillus stearothermophilus* or *Bacillus atrophaeus*. *G. stearothermophilus* has a known resistance for steam sterilization, and *B. atrophaeus* for ethylene oxide and dry heat sterilizations. When mixed with ATS or other simulated-use test soils, your reprocessing instructions from cleaning through sterilization can be validated. For more information go to:

<http://artificialtestsoil.com/bactiquik.php?pmc=PR-Bactiquik>

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.