



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Offers New EndoDolly™

Fraser, MI, August 21, 2019 - Healthmark Industries is pleased to announce the new EndoDolly™ to their Endoscope Product Line.

Manufactured from stainless steel, the EndoDolly™ is designed for hanging scopes and stands on a five-wheel locking base for ease of transportation. The EndoDolly™ is equipped with 3 self-adjusting extension poles that can be individually raised by simply touching a pressure button, which allows the user to extend the 2 outside poles to 8 feet and the middle pole to 10 feet. Located at the top of each pole is a lock tight hook holder that can accommodate two scopes per pole and supports flexible scopes up to 15 lbs. per pole. Each self-adjusting pole comes with polycarbonate cups that allow for the scope tip to be placed in during transportation, storage, or quality assurance testing.

Visit www.hmark.com or call 800-521-6224 for more information.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.