



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

New Single-Use Endoscope Cleaning Brush

Fraser, MI, June 20th 2016 -Healthmark Industries has announced the addition of the Endoscope Cleaning Brush CC-025-600 to its GI product line.

Endoscope reusable brushes are often costly and are unintentionally discarded. For this reason the single-use 6mm diameter brush has been developed as an alternate option purposed for the channels of flexible endoscopes that have a diameter between 2.0 and 3.2mm.

Endoscope Cleaning Brush is used to brush inside of the suction channel in the control section, the instrument channel junction, the forceps/irrigation plug (isolated type), the luer-split, and the inside of the suction and biopsy valves, depending on which type of endoscope is used. The nylon bristle brush tips connect to a stainless-steel wire, which joins to a versatile white thermoplastic tubing, ensuring adequate debris removal without causing scope damage.

Visit www.hmark.com or call 800-521-6224 for more information.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, support and serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information or call 800-521-6224.