



healthmark  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Offers New Instrument Inspection Tips**

Fraser, MI, August 6, 2019 - Healthmark Industries is delighted to introduce the Instrument Inspection Tips to its ProSys™ Instrument Care product line. Designed as a single-use tool, the Instrument Inspection Tips assist healthcare workers on the clean side to check if items have been effectively cleaned and are free of obstruction.

Manufactured with a spun polyester swab that connects to a flexible tubing, the Instrument Inspection Tips are offered in lengths of thirty and sixty centimeters. The Instrument Inspection Tips have the following different colored shafts that indicate the swab head diameters: White (1.1mm), Black (2.7mm) and Blue (4.2mm). The Instrument Inspection Tips are available for purchase in a package of thirty.

Visit [www.hmark.com](http://www.hmark.com) or call 800-521-6224 for more information.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.