



**healthmark**  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## **FOR IMMEDIATE RELEASE:**

### **CONTACTS:**

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Offers New Insulation Tester**

Fraser, MI May, 24th, 2019 -Healthmark Industries is delighted to introduce the Insulation Tester with optional Bi-Polar fixture to its ProSys™ Instrument Care product line.

Designed for electrically testing electrosurgical instruments, the low frequency high voltage Insulation Tester is used to detect and locate defects such as pinholes, cracks and bare spots in the jacket or coating of laparoscopic and bi-polar electrosurgical instruments. The Insulation Tester is a handheld, portable unit that tests the insulation integrity of electrosurgical equipment for flaws in protective coatings applied over conductive instrument surfaces in order to prevent inadvertent tissue burns, which may occur during electrosurgical instrument procedures.

Manufactured with a rechargeable battery, the Insulation Tester maintains applied test voltage with constant current source, features full test current at low voltages, limited output current for operational safety, easy to read LED Indicators, as well as LED display of alarm and battery charge.

The Insulation Tester comes with the following reusable and interchangeable accessories: ground wire with alligator clip, ring electrode, tri-hole electrode, brush electrode, saddle, and case. The Bi-Polar Fixture accessory that works in conjunction with the saddle, ground wire and brush electrode can be purchased separately. Additionally, an optional Wire Tester accessory used with an HV Red Lead Wire is available for testing wires by locating and identifying defects such as pinholes or cracks in the conductive core.

Visit [www.hmark.com](http://www.hmark.com) or call 800-521-6224 for more information.

### **About Healthmark Industries Co., Inc.**

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to

serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.