



**FOR IMMEDIATE RELEASE:**

CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

**Interchangeable Tags for Easy Identification**

**Fraser, MI, October 7, 2015** - Healthmark Industries has announced the addition of the Interchangeable Insert Tags to its Labeling line.

Constructed from stainless steel, the heat resistant insert holder is manufactured with a durable clip that attaches to various medical trays and mesh wire baskets. Accompanied by a variety of interchangeable color insert tags for easy identification, these can be laser imprinted with text and barcode (green labels can only be engraved with text). Made of resilient polyethersulfone, the interchangeable insert tags are heat resistant up to 275 °F.

**About Healthmark Industries Co., Inc.**

*Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, support and serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information or call 800-521-6224.*