



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Offers new LED Interchangeable Magnifier

Fraser, MI, November 19, 2019 - Healthmark Industries has announced the addition of the LED Interchangeable Magnifier to its ProSys™ Instrument Care product line.

Designed for items to be optically inspected, the LED Interchangeable Magnifier is manufactured with an interchangeable 5-diopter glass lens that has a 2.25x magnification, 360 rotatable lamp head and a fully enclosed spring balanced arm to prevent pinching hazards. Equipped with 64 ultra-bright white SMD LEDs with brightness adjustable controls, the LED Interchangeable Magnifier has a reach of 36 Inches, max illumination of 950lm, wattage of 6W and includes a table clamp with a 2-inch opening. The LED Interchangeable Magnifier power supply has an input of 100-240V-50/60Hz -0.45A and output of 26V/500mA.

Interchangeable lenses are offered in the following magnifications and diopters: 2.25x (*5 diopter*), 3x (*8 diopter*), 4x (*12 diopter*), and 4.75x (*15 diopter*). The LED Interchangeable Magnifier is available for individual purchase. The Lenses can be purchased separately or together

Visit www.hmark.com or call 800-521-6224 for more information.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.