



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Offers New Cable Continuity Tester

Fraser, MI October 16, 2020 - Healthmark Industries is delighted to introduce the Cable Continuity Tester to its ProSys™ Instrument Care product line.

Designed to test the quality of monopolar and bipolar cords, the Cable Continuity Tester is user-friendly and portable. For the desired test, simply position the “ON” switch to bipolar or monopolar. If the green light remains illuminated, then the cable passed testing. If the light flickers or doesn’t illuminate, this indicates the cable failed testing. Each Cable Continuity Tester comes with a 9V battery and a carrying case.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries’ mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.