

## Instructions for Use: Beard Net

<b>Brand Name of Product</b>	Beard Net
<b>Generic Name of Product</b>	Beard Net
<b>Product Code Number(s)</b>	BN300
<b>Intended Use</b>	To cover the beard to limit as a source of possible contamination in the prep and pack area and decontamination area.
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>● One size fits most</li> <li>● Super fine mesh nylon with small holes</li> <li>● Single-use</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	100 in a box
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	To cover the beard to limit as a source of possible contamination.
<b>Preparation for Use</b>	Take one beard net from the box.
<b>Diagrams (drawings, pictures)</b>	N/A
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Place hands around the band of the beard net.</li> <li>2. Place the beard net over the head and around facial hair.</li> <li>3. The band is placed around the back of the head and above the ears.</li> <li>4. Make sure all facial hair is completely covered by the beard net.</li> <li>5. To remove, with both hands take the band and lift over the head and pull the beard net away from face.</li> </ol>
<b>Interpretation of Test Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>● One person, one use</li> <li>● Not intended for use in the procedure areas.</li> </ul>
<b>Disposal</b>	Single-use only

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	PPE Brochure, PPE Price List

<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2021-04-23 Suzanne Latta