

Instructions for Use: Container Loks without Indicator Dot



<b>Brand Name of Product</b>	Container Loks
<b>Generic Name of Product</b>	Closed sterilization container tamper evident seal
<b>Product Code Number(s)</b>	3030 BL, 3030 WT, 3030 YL, 3030 OR, 4040 GN, 605LNG GR, 605M GR
<b>Intended Use</b>	Use the tamper-evident seal to indicate the container has not been opened since sterilization
<b>Range of Applications for Product</b>	N/A
<b>Key Specifications of Product</b>	Blue, Yellow, White, Orange, Green (4040 only available in green)

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

**Instructions for Using Product**

<b>Description of Use(s)</b>	Use the tamper-evident seal to indicate the container has not been opened since sterilization.
------------------------------	--

<b>Preparation for Use</b>	N/A
----------------------------	-----

<b>Diagrams (drawings, pictures)</b>	   <p style="text-align: center;"><b>Container Lok without Indicator Dot</b></p>
--------------------------------------	--

<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Pass arrow shackle through opening of item to be sealed.</li> <li>2. Hold the base of the lok with one hand and firmly push down on the arrow shackle into locking device with the other hand until it engages/locks with a small snap sound. Use one hand to hold the base of the lok and the other hand to pull on seal to ensure it has been properly locked.</li> </ol>
---------------------------------	---

<b>Interpretation of Test Results</b>	N/A
---------------------------------------	-----

<b>Contraindications of Test Results</b>	N/A
--	-----

<b>Documentation</b>	N/A
----------------------	-----

<b>Special Warnings and Cautions</b>	If seal is broken, assume sterility has been compromised and the tray must be reprocessed.
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	Tamper-evident seal is materially compatible with steam, VH2O2/O3/EO sterilization.
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	Sterilization Packaging Products Brochure, Sterilization Packaging Brochure Products Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com