

Instructions for Use: Cool Aids Single-Use Vest

Brand Name of Product	Cool Aids Single-Use Vest
Generic Name of Product	Cooling Vest
Product Code Number(s)	SUCV-001, SUCV-002
Intended Use	To provide cooling of the core body temperature.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> One person, one vest per use Includes inside pockets for phase change packs

Shipping & Storage

Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product

Description of Use (s)	To provide cooling of the core body temperature.
Preparation for Use	Activate by placing phase change inserts in a refrigerator or freezer to cool. It is recommended to wear a thin undergarment such as a cotton T-Shirt under the vest for comfort. Simply place CoolPax™ inserts in freezer for 35/40 minutes and insert into garment.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> Place each phase change packs inside the mesh pockets inside the vest. Wait 10 minutes before putting on the vest. Slide vest overhead so that it rests on shoulders. Tie side straps on each side.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	It is recommended to dispose of in compliance with the facility's protocols for worn PPE.

Reprocessing Instructions

Point of use	N/A
Preparation for decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	<u>Cooling inserts</u> <ul style="list-style-type: none"> Clean with light soap and water
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, inspection, and testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	PPE Brochure, PPE Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com

