



Instructions for Use: Dual UnderGuard™ Dry Mats

Brand Name of Product	Dual UnderGuard™ Dry Mats
Generic Name of Product	Dry Mats
Product Code Number(s)	DM-2030PP, DM-2030GN, DM-1520PP, DM-1520GN, DM-1224PP, DM-1224GN, DM-1020PP, DM-1220PP, DM-0820PP, DM-1212GN, DM-1212PP
Purpose of Product	To be used under trays and inside of trays with the purple side facing down to protect against handling while simultaneously providing an absorbent inner layer for the dispersion of condensate.
Range of Applications for Product	Provides protection useful for heavy trays and baskets.
Key Specifications of Product	Dry mats are made from polypropylene and viscose.

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	+15°C-30°C and 30-70% RH
Packaging Contents	
Shelf Life	

Instructions for Using Product	
Description of Use(s)	Strong spun polypropylene sheet with an absorbent layer to protect against sharp-edged trays and promote the dispersion of condensate.
Preparation for Use	
Diagrams (drawings, pictures)	
Steps for Use of Product	<ul style="list-style-type: none"> Liner may be placed either under the tray, or inside the bottom of the tray.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	Dry mats are not compatible with VH2O2 sterilization processes.
Disposal	Single-use. Dispose of according to facility's policy.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	<ul style="list-style-type: none"> Steam sterilization at 132°C for 4 minutes or 135°C for 3 minutes Gravity steam sterilization at 121°C for 30 minutes Maximum tested time/temperature: 30 minutes at 135°C Materially compatible with EO sterilization
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	Sterilization Brochure, Sterilization Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. Fraser, MI 48026 1-586-774-7600

healthmark@hmark.com
hmark.com