

Instructions for Use: HTK Labels

Brand Name of Product	HangTimeKit™ Labels
Generic Name of Product	Kit for managing and insuring the safe storage of flexible endoscopes
Product Code Number(s)	403225 HTK, 403225 HTKG, 403225 HTKY, 403225 HTKB, 403225 HTKO, 403225 HTKP
Purpose of Product	The HangTime label is for documenting the date of reprocessing.
Range of Applications for Product	Surveillance of flexible endoscopes in storage to track the date of when the endoscope was last reprocessed.
Key Specifications of Product	<ul style="list-style-type: none"> Numbers representing the days of the month and abbreviations for the 12 months are printed on the label. Offered in a variety of colors.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Conditions	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use (s)	The label is for recording and documenting the latest date that the flexible endoscope was reprocessed.
Preparation for Use	N/A
Diagrams (drawings, pictures):	N/A
Steps for Use of Product	<ol style="list-style-type: none"> Circle the day and month of the last reprocessed date. Optionally include the initials of the reprocessor. Place the label securely around the scope through the looped opening.
Interpretation of Results	If the endoscope has been in storage for an extended amount of time, you may need to test it with the HangTimeKit™ to be sure it's still safe to use.
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	Labels are single use only. When finished, it should be treated as biohazard and disposed of in accordance to your facilities guidelines.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	HangTimeKit™
Other Product Support Documents	Hang Time Kit Specification Sheet, HangTimeKit Pocket Brochure, HangTimeKit Whitepaper, HTK-FV1041 Gram Negative Test Specification Sheet, Illustrated Users Guide for HangTime test
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc

	18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com
--	--