

Brand Name of Product	Nylon Cart Covers
Generic Name of Product	Re-useable Cart Covers
Product Code Number(s)	SRC-835, SRC-836, SRC-837, SRC-435, SRC-436, SRC-437, SRC-438, SRC-439, SRC-440, SRC-441, SRC-442, SRC-443, SRC-444
Intended Use	Ideal for transport carts with an extra level of protection.
Range of Applications for Product	Medical supplies and accessories
Key Specifications of Product	200 Denier Nylon: <ul style="list-style-type: none"> • Woven material for stationary carts • Machine washable • Uncoated covers 400 Denier Nylon: <ul style="list-style-type: none"> • Tight woven for transport carts (i.e.: linens and laundry carts) • Machine washable • Uncoated covers

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Ideal for transport carts with an extra level of protection.
Preparation for Use	N/A
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	Place cover over the cart so it is completely covered.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	<ul style="list-style-type: none"> • Do not put covers through flatwork ironers • Washing or drying more than recommended in the <i>Cleaning-Manual</i> section or <i>Drying</i> section will reduce service life and performance of the product
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	<ul style="list-style-type: none"> • Use mild non-bleach detergent suitable for nylon/polyester fabrics and follow manufacturer's instructions • Water temperature not to exceed 71 °C (160 °F)
Disinfection	N/A
Drying	Air dry. If machine drying is necessary, do not exceed a temperature of 132 °C (250 °F) maximum and time duration of 2-minutes maximum.
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A

Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	Bins Galore Brochure, Bins Galore Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com