

	A	B	C	D	E	F	G	H
	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
1								
2	Laparoscopic Trays							
3		DDL-1123	Double Decker Lap Tray	Anodized Aluminum	Y	Y	Y	
4		CYSTO-2214B	Cysto Tray	Stainless Steel	Y	Y	Y	
5		DDL-4113	Poly Double Decker Tray	Radel	Y	Y	Y	
6		51527	Double Decker Bariatric	Anodized Aluminum	Y	Y	Y	
7		31127 DV	Da Vinci Arm Sterilization Case	Anodized Aluminum	Y	Y	Y	
8		300605-DVC	Da Vinci Lens Tray	Anodized Aluminum	Y	Y	Y	
9		300605-DVA	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
10		300605-H	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
11		180502-INS	Da Vinci Tray Insert Arm Case	Anodized Aluminum	Y	Y	Y	
12								
13	Custom Procedure Trays							
14		30320-H	2 Level Tray System	Anodized Aluminum	Y	Y	Y	
15		31421-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
16		31123-T	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
17		31123-T1	Universal Tray-Double Layer	Anodized Aluminum	Y	Y	Y	
18		31127-H	Universal Tray-Single Layer	Anodized Aluminum	Y	Y	Y	
19								
20	Instrument Brackets for Custom Trays							
21		P-01	Post	Anodized Aluminum/Silicone	Y	Y	Y	
22		Tri-CUT BRACKET	1"universal bracket 3, 5 & 10 mm	Stainless Steel/Silicone	Y	Y	Y	
23		56-UB	5" 6 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
24		710-UB	7.5" 10 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
25		1015-UB	10" 15 position universal bracket	Anodized Aluminum/Silicone	Y	Y	Y	
26								
27	Instru-Trays							
28		4116-H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
29		4106H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
30		41013H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
31		41116H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
32		41317H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
33		41520H	Perforated Tray	Polypropylene	Y	Y	Y	May warp slightly after cleaning and/or sterilization, performance not affected.
34		HCO-1309	Perforated Tray	Stainless Steel	Y	Y	Y	
35		HCO-1510	Perforated Tray	Stainless Steel	Y	Y	Y	
36		HCO-1711	Perforated Tray	Stainless Steel	Y	Y	Y	
37		HCO-1912	Perforated Tray	Stainless Steel	Y	Y	Y	
38		HCO-1065	Perforated Tray	Stainless Steel	Y	Y	Y	
39		HCO-2116	Perforated Tray	Stainless Steel	Y	Y	Y	
40		HCO-1065-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
41		HCO-1309-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
42		HCO-1510-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
43		HCO-1711-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
44		HCO-1912-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
45								
46	Metal Micro Trays							
47		184014	Base & Cover	Anodized Aluminum	Y	Y	Y	
48		184016	Base & Cover	Anodized Aluminum	Y	Y	Y	
49		184019	Base & Cover	Anodized Aluminum	Y	Y	Y	
50		184021	Base & Cover	Anodized Aluminum	Y	Y	Y	
51		181453	Solid (no holes) Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
52		184453	Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
53		181463	Base,Cover & Holder	Anodized Aluminum	Y	Y	Y	
54		182500	Stainless Steel - 5 Clip	Stainless Steel	Y	Y	Y	
55		182062	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
56		182068	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
57								
58	Hinged Top Micro Trays							
59		30732-H	Taupe	Polypropylene	Y	Y	Y	
60		30842-H	Taupe	Polypropylene	Y	Y	Y	
61		31042-H	Taupe	Polypropylene	Y	Y	Y	
62		31055-H	White	Polypropylene	Y	Y	Y	
63		31462-H	Taupe	Polypropylene	Y	Y	Y	
64		31772-H	Blue	Polypropylene	Y	Y	Y	
65		50Z900	White	Polypropylene	Y	Y	Y	
66								
67	Hinged Top Non-Perforated Trays							
68		30732	Taupe	Propylene/Ethylene	Y	Y	Y	
69		30842	Taupe	Propylene/Ethylene	Y	Y	Y	
70		31042	Taupe	Propylene/Ethylene	Y	Y	Y	
71		31462	Taupe	Propylene/Ethylene	Y	Y	Y	
72		31772	Blue	Propylene/Ethylene	Y	Y	Y	
73								
74	Specialty Micro Trays							
75		30109HN	ER Hand Trays	Anodized Aluminum	Y	Y	Y	
76		30109CC	Camera Case	Anodized Aluminum	Y	Y	Y	
77		30011KE	Key Elevators Box	Anodized Aluminum	Y	Y	Y	
78		30011CU	Currettes Box	Anodized Aluminum	Y	Y	Y	
79		30011RS	Rhoton Suction Box	Stainless Steel	Y	Y	Y	
80		30011BS	Bipolar Set Box	Stainless Steel	Y	Y	Y	
81		30011BP	Bipolar Set Box	Stainless Steel	Y	Y	Y	
82								
83	Metal Hinged Top Micro Trays							
84		31010-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
85		31010-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
86		31015-H1	Tray & Cover	Anodized Aluminum	Y	Y	Y	
87		31015-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
88		31020-H	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
89								
90	Clear-Top Micro Trays							
91		0278	Base & Cover	Ultem	Y	Y	Y	
92		02781	Base & Cover	Ultem	Y	Y	Y	
93		027814	Base & Cover	Ultem	Y	Y	Y	
94		0625	Base & Cover	Ultem	Y	Y	Y	
95		062514	Base & Cover	Ultem	Y	Y	Y	
96		6540	Base & Cover	Ultem	Y	Y	Y	
97		7540	Base & Cover	Ultem	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
98		754015	Base & Cover	Ultem	Y	Y	Y	
99		100675	Single Decker Twin	Ultem	Y	Y	Y	
100		100615	Double Decker Twin	Ultem	Y	Y	Y	
101		100615B	Insert Tray for 100615	Ultem	Y	Y	Y	
102		101525	Base & Lid	Ultem	Y	Y	Y	
103		101525B	Full Insert	Ultem	Y	Y	Y	
104								
105	Silicone Plug-InSecure-It							
106		A237591K	One 4mm slot, two plugs	Silicone	Y	Y	Y	
107		A237592K	Two 4mm slots, two plugs	Silicone	Y	Y	Y	
108		A237593K	Three 4mm slots, two plugs	Silicone	Y	Y	Y	
109		A237593	Three 1mm slots, two plugs	Silicone	Y	Y	Y	
110		A237596	Six 1mm slots, two plugs	Silicone	Y	Y	Y	
111		A387596	Six 1mm slots, three plugs	Silicone	Y	Y	Y	
112		A3875912	Twelve 1mm slots, three plugs	Silicone	Y	Y	Y	
113		A387592	Two 4mm slots, three plugs	Silicone	Y	Y	Y	
114		A387593	Three 4mm slots,three plugs	Silicone	Y	Y	Y	
115		A387594	Four 4mm slots, three plugs	Silicone	Y	Y	Y	
116		A387595	Five 4mm slots, three plugs	Silicone	Y	Y	Y	
117		A2616901	One 12mm slots, two plugs	Silicone	Y	Y	Y	
118		A26216902	Two 8mm slots, two plugs	Silicone	Y	Y	Y	
119		A32759P2	One 5mm slot, three plugs	Silicone	Y	Y	Y	
120		A327593	Three 5mm slots, three plugs	Silicone	Y	Y	Y	
121		A327593L	Two 12mm slots, three plugs	Silicone	Y	Y	Y	
122		A32759P1	One 12mm slots, three plugs	Silicone	Y	Y	Y	
123								
124	Scope Trays							
125		30126	Bariatric - Single Scope Case	Anodized Aluminum	Y	Y	Y	
126		30630	Bariatric - Multi Scope Case	Anodized Aluminum	Y	Y	Y	
127		30111-SSC	Scope Tray	Stainless Steel	Y	Y	Y	
128		30118-SSC	Universal Scope Case	Stainless Steel	Y	Y	Y	
129		31042-22	Hinged Top Arthro Tray	Polypropylene	Y	Y	Y	
130		3-4110	Clear Top Arthro Tray	Radel	Y	Y	Y	
131		SP-3014	Osteotome Tray	Anodized Aluminum	Y	Y	Y	
132		3-5120	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
133		3-5130	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
134		75.0047.9M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
135		75.0083.1M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
136		75.0047.5M	Wire Basket Scope Tray	Stainless Steel	Y	Y	Y	
137		750083	Wire Scope Tray	Stainless Steel	Y	Y	Y	
138		31772-22	Economy Scope Tray	Polypropylene	Y	Y	Y	
139		261615	Clear Top Scope Tray	Radel	Y	Y	Y	
140		302215	Scope Tray B/L 1 set Secur-Its	Radel	Y	Y	Y	
141		31462-22	Camera Tray	Polypropylene	Y	Y	Y	
142		3-3131	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
143		3-3170	One Scope 1mm-10mm	Radel	Y	Y	Y	
144		3-3180	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
145		3-3200	One Scope 1mm-10mm	Radel	Y	Y	Y	
146		9326	Super Long Scope Tray	Radel	Y	Y	Y	
147		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
148								
149								
150	Silicone Finger Mats							
151		202025	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
152		2-1405	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
153		SP2313	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
154		SP4622	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
155		SP2719	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
156		2-3215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
157		2-1055	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
158		2-1215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
159		6-1914	Silicone Mat Bubble Bottom (no holes)	Silicone	Y	Y	Y	
160		2-4005	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
161		2111-M	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
162		2-6005	Silicone Mat-Bubble Bottom (2 piece)	Silicone	Y	Y	Y	
163		SR10	Silicone Mat-Perforated Base - Roll -Pink	Silicone	Y	Y	Y	
164		0278-M	Silicone Mat 2.667 x 1.5	Silicone	Y	Y	Y	
165		0625-M	Silicone Mat 2.5 W x 6L	Silicone	Y	Y	Y	
166		1006158M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
167		100615-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
168		100675-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
169		101525B-M	Finger Mat	Silicone	Y	Y	Y	
170		6540-M	Micro-Protech Mat 4 x 7.5	Silicone	Y	Y	Y	
171		7540-M	Silicone Mat	Silicone	Y	Y	Y	
172								
173	Wire Baskets							
174		75.0001.7	Wire Basket	Stainless Steel	Y	Y	Y	
175		75.0005.3	Wire Basket	Stainless Steel	Y	Y	Y	
176		75.0004.6	Wire Basket	Stainless Steel	Y	Y	Y	
177		75.0005.9	Wire Basket	Stainless Steel	Y	Y	Y	
178		75.0004.8	Wire Basket	Stainless Steel	Y	Y	Y	
179		75.0004.9	Wire Basket	Stainless Steel	Y	Y	Y	
180		75.0072.9	Wire Basket	Stainless Steel	Y	Y	Y	
181		75.0091.1	Cover	Stainless Steel	Y	Y	Y	
182		75.0091.2	Cover	Stainless Steel	Y	Y	Y	
183		7501133	Wire Basket with Lid	Stainless Steel	Y	Y	Y	
184		72.1400.5	Large Washer Basket	Stainless Steel	Y	Y	Y	
185		72.1400.6	Large Washer Basket	Stainless Steel	Y	Y	Y	
186								
187	High-Sided Perforated Trays							
188		HDO-0952	Perforated Tray	Stainless Steel	Y	Y	Y	
189		HDO-1062	Perforated Tray	Stainless Steel	Y	Y	Y	
190		HDO-1272	Perforated Tray	Stainless Steel	Y	Y	Y	
191		HDO-1210	Perforated Tray	Stainless Steel	Y	Y	Y	
192		HDO-1610	Perforated Tray	Stainless Steel	Y	Y	Y	
193		HDO-2012	Perforated Tray	Stainless Steel	Y	Y	Y	
194								
195	High-Sided Non-Perforated Trays							
196		HDO-1062-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
197		HDO-1210-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
198		HDO-1272-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
199		HDO-1610-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
200		HDO-2012-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
201		HDO-952-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
202								
203	Plastic Sterilization Baskets							
204		PIT2015	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
205		DINT2525	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
206		4106-H	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
207		PIT2715	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
208		MIT30251	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
209		1511-4	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
210		DT-1218	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
211		MIT40301	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
212		DT-2136	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
213		1511-7	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
214		DINT4825	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
215		2006-H	Plastic Sterilization Baskets	Polyarylethersulfone	Y	Y	Y	
216		113-1314H	Plastic Sterilization Baskets	Polyarylethersulfone	Y	Y	Y	
217								
218	Perforated Stainless Steel Trays							
219		PSST-009	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
220		PSST-007	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
221		PSST-005	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
222		PSST-010	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
223		PSST-001	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
224		PSST-002	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
225		PSST-012	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
226		PSST-003	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
227		PSST-011	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
228		PSST-004	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
229		PSST-013	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
230		PSST-008	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
231		PSST-006	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
232								
233	Long Scope Instrument Tray							
234		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
235								
236	Sterilizing Basket							
237		PIT3025	Sterilizing Basket 11.82 x 9.85 x 2.05H	Stainless Steel	Y	Y	Y	
238		PIT5739	Sterilizing Basket 21.81 x 15.44 x 2.05H	Stainless Steel	Y	Y	Y	
239								
240	Nic-Nac Washer Baskets							
241		CIL-004	Fine Mesh Basket with Snap Closure	Stainless Steel	Y	Y	Y	
242		CIL-005	Fine Mesh Basket with Snap Closure	Stainless Steel	Y	Y	Y	