



	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH202/O3 Compatibility	Additional Comments
45	Metal Micro Trays							
46		184014	Base & Cover	Anodized Aluminum	Y	Y	Y	
47		184016	Base & Cover	Anodized Aluminum	Y	Y	Y	
48		184019	Base & Cover	Anodized Aluminum	Y	Y	Y	
49		184021	Base & Cover	Anodized Aluminum	Y	Y	Y	
50		181453	Solid (no holes) Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
51		184453	Tray w/ Holder	Anodized Aluminum	Y	Y	Y	
52		181463	Base,Cover & Holder	Anodized Aluminum	Y	Y	Y	
53		182500	Stainless Steel - 5 Clip	Stainless Steel	Y	Y	Y	
54		182062	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
55		182068	Stainless Steel Cassette	Stainless Steel	Y	Y	Y	
56								
57	Hinged Top Micro Trays							
58		30732-H	Taupe	Polypropylene	Y	Y	Y	
59		30842-H	Taupe	Polypropylene	Y	Y	Y	
60		31042-H	Taupe	Polypropylene	Y	Y	Y	
61		31055-H	White	Polypropylene	Y	Y	Y	
62		31462-H	Taupe	Polypropylene	Y	Y	Y	
63		31772-H	Blue	Polypropylene	Y	Y	Y	
64		50Z900	White	Polypropylene	Y	Y	Y	
65								
66	Hinged Top Non-Perforated Trays							
67		30732	Taupe	Propylene/Ethylene	Y	Y	Y	
68		30842	Taupe	Propylene/Ethylene	Y	Y	Y	
69		31042	Taupe	Propylene/Ethylene	Y	Y	Y	
70		31462	Taupe	Propylene/Ethylene	Y	Y	Y	
71		31772	Blue	Propylene/Ethylene	Y	Y	Y	
72								
73	Specialty Micro Trays							
74		30109HN	ER Hand Trays	Anodized Aluminum	Y	Y	Y	
75		30109CC	Camera Case	Anodized Aluminum	Y	Y	Y	
76		30011KE	Key Elevators Box	Anodized Aluminum	Y	Y	Y	
77		30011CU	Currettes Box	Anodized Aluminum	Y	Y	Y	
78		30011RS	Rhoton Suction Box	Stainless Steel	Y	Y	Y	
79		30011BS	Bipolar Set Box	Stainless Steel	Y	Y	Y	
80		30011BP	Bipolar Set Box	Stainless Steel	Y	Y	Y	
81								
82	Metal Hinged Top Micro Trays							
83		31015-H2	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
84		31020-H	Tray & Cover & Insert	Anodized Aluminum	Y	Y	Y	
85								
86	Clear-Top Micro Trays							
87		0278	Base & Cover	Ultem	Y	Y	Y	
88		02781	Base & Cover	Ultem	Y	Y	Y	
89		027814	Base & Cover	Ultem	Y	Y	Y	
90		0625	Base & Cover	Ultem	Y	Y	Y	
91		062514	Base & Cover	Ultem	Y	Y	Y	
92		6540	Base & Cover	Ultem	Y	Y	Y	
93		7540	Base & Cover	Ultem	Y	Y	Y	
94		754015	Base & Cover	Ultem	Y	Y	Y	
95		100675	Single Decker Twin	Ultem	Y	Y	Y	
96		100615	Double Decker Twin	Ultem	Y	Y	Y	
97		100615B	Insert Tray for 100615	Ultem	Y	Y	Y	
98		101525	Base & Lid	Ultem	Y	Y	Y	
99		101525B	Full Insert	Ultem	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
100								
101	Silicone Plug-InSecure-It							
102		A237591K	One 4mm slot, two plugs	Silicone	Y	Y	Y	
103		A237592K	Two 4mm slots, two plugs	Silicone	Y	Y	Y	
104		A237593K	Three 4mm slots, two plugs	Silicone	Y	Y	Y	
105		A237593	Three 1mm slots, two plugs	Silicone	Y	Y	Y	
106		A237596	Six 1mm slots, two plugs	Silicone	Y	Y	Y	
107		A387596	Six 1mm slots, three plugs	Silicone	Y	Y	Y	
108		A3875912	Twelve 1mm slots, three plugs	Silicone	Y	Y	Y	
109		A387592	Two 4mm slots, three plugs	Silicone	Y	Y	Y	
110		A387593	Three 4mm slots,three plugs	Silicone	Y	Y	Y	
111		A387594	Four 4mm slots, three plugs	Silicone	Y	Y	Y	
112		A387595	Five 4mm slots, three plugs	Silicone	Y	Y	Y	
113		A2616901	One 12mm slots, two plugs	Silicone	Y	Y	Y	
114		A26216902	Two 8mm slots, two plugs	Silicone	Y	Y	Y	
115		A32759P2	One 5mm slot, three plugs	Silicone	Y	Y	Y	
116		A327593	Three 5mm slots, three plugs	Silicone	Y	Y	Y	
117		A327593L	Two 12mm slots, three plugs	Silicone	Y	Y	Y	
118		A32759P1	One 12mm slots, three plugs	Silicone	Y	Y	Y	
119								
120	Scope Trays							
121		30126	Bariatric - Single Scope Case	Anodized Aluminum	Y	Y	Y	
122		30630	Bariatric - Multi Scope Case	Anodized Aluminum	Y	Y	Y	
123		30111-SSC	Scope Tray	Stainless Steel	Y	Y	Y	
124		30118-SSC	Universal Scope Case	Stainless Steel	Y	Y	Y	
125		31042-22	Hinged Top Arthro Tray	Polypropylene	Y	Y	Y	
126		3-4110	Clear Top Arthro Tray	Radel	Y	Y	Y	
127		SP-3014	Osteotome Tray	Anodized Aluminum	Y	Y	Y	
128		3-5120	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
129		3-5130	Scope Tray w/Cord Channel	Radel	Y	Y	Y	
130		750083	Wire Scope Tray	Stainless Steel	Y	Y	Y	
131		31772-22	Economy Scope Tray	Polypropylene	Y	Y	Y	
132		261615	Clear Top Scope Tray	Radel	Y	Y	Y	
133		302215	Scope Tray B/L 1 set Secur-Its	Radel	Y	Y	Y	
134		31462-22	Camera Tray	Polypropylene	Y	Y	Y	
135		3-3131	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
136		3-3170	One Scope 1mm-10mm	Radel	Y	Y	Y	
137		3-3180	Two Scopes 1mm-10mm	Radel	Y	Y	Y	
138		3-3200	One Scope 1mm-10mm	Radel	Y	Y	Y	
139		9326	Super Long Scope Tray	Radel	Y	Y	Y	
140		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
141								
142								
143	Silicone Finger Mats							
144		202025	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
145		2-1405	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
146		SP2313	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
147		SP4622	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
148		SP2719	Silicone Mat-Perforated Base	Silicone	Y	Y	Y	
149		2-3215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
150		2-1055	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
151		2-1215	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
152		6-1914	Silicone Mat Bubble Bottom (no holes)	Silicone	Y	Y	Y	
153		2-4005	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	
154		2111-M	Silicone Mat Bubble Bottom	Silicone	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH2O2/O3 Compatibility	Additional Comments
155		2-6005	Silicone Mat-Bubble Bottom (2 piece)	Silicone	Y	Y	Y	
156		SR10	Silicone Mat-Perforated Base - Roll -Pink	Silicone	Y	Y	Y	
157		0278-M	Silicone Mat 2.667 x 1.5	Silicone	Y	Y	Y	
158		0625-M	Silicone Mat 2.5 W x 6L	Silicone	Y	Y	Y	
159		1006158M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
160		100615-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
161		100675-M	Silicone Mat 6W x 10L	Silicone	Y	Y	Y	
162		101525B-M	Finger Mat	Silicone	Y	Y	Y	
163		6540-M	Micro-Protech Mat 4 x 7.5	Silicone	Y	Y	Y	
164		7540-M	Silicone Mat	Silicone	Y	Y	Y	
165								
166	Wire Baskets							
167		75.0005.9	Wire Basket	Stainless Steel	Y	Y	Y	
168		75.0072.9	Wire Basket	Stainless Steel	Y	Y	Y	
169		75.0091.2	Cover	Stainless Steel	Y	Y	Y	
170		72.1400.5	Large Washer Basket	Stainless Steel	Y	Y	Y	
171		72.1400.6	Large Washer Basket	Stainless Steel	Y	Y	Y	
172								
173	High-Sided Perforated Trays							
174		HDO-0952	Perforated Tray	Stainless Steel	Y	Y	Y	
175		HDO-1062	Perforated Tray	Stainless Steel	Y	Y	Y	
176		HDO-1272	Perforated Tray	Stainless Steel	Y	Y	Y	
177		HDO-1210	Perforated Tray	Stainless Steel	Y	Y	Y	
178		HDO-1610	Perforated Tray	Stainless Steel	Y	Y	Y	
179		HDO-2012	Perforated Tray	Stainless Steel	Y	Y	Y	
180								
181	High-Sided Non-Perforated Trays							
182		HDO-1062-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
183		HDO-1210-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
184		HDO-1272-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
185		HDO-1610-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
186		HDO-2012-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
187		HDO-952-NP	Non-Perforated Tray	Stainless Steel	Y	Y	Y	
188								
189	Plastic Sterilization Baskets							
190		PIT2015	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
191		DINT2525	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
192		4106-H	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
193		PIT2715	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
194		MIT30251	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
195		1511-4	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
196		DT-1218	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
197		MIT40301	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
198		DT-2136	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
199		1511-7	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
200		DINT4825	Plastic Sterilization Baskets	Polypropylene	Y	Y	Y	
201		2006-H	Plastic Sterilization Baskets	Polyetherimide	Y	Y	Y	
202		113-1314H	Plastic Sterilization Baskets	Polyetherimide	Y	Y	Y	
203								
204	Perforated Stainless Steel Trays							
205		PSST-009	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
206		PSST-007	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
207		PSST-005	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
208		PSST-010	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
209		PSST-001	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	

	A	B	C	D	E	F	G	H
1	System Code No	Component Parts	Description	Material Composition	Steam Sterilization Compatibility	Ethylene Oxide Compatibility	VH202/03 Compatibility	Additional Comments
210		PSST-002	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
211		PSST-012	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
212		PSST-003	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
213		PSST-011	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
214		PSST-004	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
215		PSST-013	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
216		PSST-008	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
217		PSST-006	Perforated Stainless Steel Tray	Stainless Steel	Y	Y	Y	
218								
219	Long Scope Instrument Tray							
220		SMT5510	Scope Mesh Container	Polypropylene	Y	Y	Y	
221								
222	Sterilizing Basket							
223		PIT3025	Sterilizing Basket 11.82 x 9.85 x 2.05H	Stainless Steel	Y	Y	Y	
224		PIT5739	Sterilizing Basket 21.81 x 15.44 x 2.05H	Stainless Steel	Y	Y	Y	
225								
226	Nic-Nac Washer Baskets							
227		CIL-004	Fine Mesh Basket with Snap Closure	Stainless Steel	Y	Y	Y	
228		CIL-005	Fine Mesh Basket with Snap Closure	Stainless Steel	Y	Y	Y	