


Instructions for Use: Silicone Tray Corner Protector

Brand Name of Product	Silicone Tray Corner Protectors
Generic Name of Product	Tray Corner Protector
Product Code Number(s)	CP-80NS, SCH-01, SCH-02
Purpose of Product	To protect the corner of the sterilization wrap from the sharp edges of sterilization trays.
Range of Applications for Product	<ul style="list-style-type: none"> For wrapped trays that are steam, EO, or gas plasma sterilized. Tray corners are reusable.
Key Specifications of Product	<ul style="list-style-type: none"> Material made of silicone. Colors (per product code): <ul style="list-style-type: none"> CP-80NS: Pink SCH-01: Purple SCH-02: Yellow

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	These reusable tray corners are a perfect solution for problematic corner tears.
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	Place tray on the tray corner protectors. Ensure all four corners are covered.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Tray Corner Protector can be cleaned manually with detergent and water.
Cleaning – Automated	Tray Corner Protectors can be processed through an automated instrument washer.
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A

Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> Compatible with steam sterilization: <ul style="list-style-type: none"> 132 °C (270 °F) for four (4) minutes 135 °C (275 °F) for three (3) minutes. VH2O2/O3/EO sterilization.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Sterilization Brochure, Sterilization Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com