



<b>Brand Name of Product</b>	Snap-Loks
<b>Generic Name of Product</b>	Tamper Evident Seals
<b>Product Code Number(s)</b>	7525, 8050
<b>Purpose of Product</b>	To show tampering and/or access to critical equipment and supplies.
<b>Range of Applications for Product</b>	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Numbered and Unnumbered</li> <li>• Seal and broken by hand</li> <li>• Cannot be resealed once broke</li> <li>• Variety of colors available</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	100 per box
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use (s)</b>	Seal crash carts, emergency boxes and other areas where evidence of tampering is critical
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	<p>The diagrams illustrate the following steps:</p> <ol style="list-style-type: none"> <li>1. Inserting the seal into the opening of the object.</li> <li>2. Pushing the seal through the locking dome until it is seated.</li> <li>3. Pressing the two sides of the object together to snap the seal shut.</li> <li>4. The final sealed state with the seal bridging the two sides.</li> </ol>
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Insert seal into closing of two sides of the object being secured.</li> <li>2. Seal through the locking dome of the seal. Snap item together until Lok is tightly sealed.</li> <li>3. The seal can be easily opened by breaking by hand by applying pressure at the scored point on the hasp.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	

<b>Documentation</b>	Follow the policy of your facility for recording date, staff and other information critical to management of critical inventory.
<b>Special Warnings and Cautions</b>	If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory.
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	Compatible with steam sterilization. Standard cycle: 270 °F/4 minutes
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Security Products Brochure, Security Products Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>