



Instructions for Use: Transport Trolley

Brand Name of Product	Transport Trolley
Generic Name of Product	Transport Trolley
Product Code Number(s)	2220, 5-2220
Intended Use	To retrieve round soaker trays and transport them for reprocessing.
Range of Applications for Product	N/A
Key Specifications of Product	N/A

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To retrieve round soaker trays and transport them for reprocessing.

Preparation for Use

Parts Included for Assembly

Figure 1

Diagrams (drawings, pictures)

	Figure 2	Figure 3	Figure 4
			
	Figure 5	Figure 6	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Screw bolt (I), according to the <i>Parts Included for Assembly</i> drawing (Fig. 1), using tool wrench (M) and place the cover cap (K). (Fig. 2). 2. Assemble frame as shown in (Fig. 3), using a plastic hammer to prevent frame from getting damaged (C). Assemble top shelf (R) using wrench (L) to screw the bolts (H) and washer (J) onto the trolley. (Fig. 4). 		
Interpretation of Results	N/A		
Contraindications of Test Results	N/A		
Documentation	N/A		
Special Warnings and Cautions	N/A		
Disposal	N/A		

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Wipe down trolley with a cleaning detergent using a moist towel.
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	PPE Brochure, PPE Price List
Other Product Support Documents	N/A
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com