


Instructions for Use: Twist-Loks

Brand Name of Product	Twist-Loks
Generic Name of Product	Tamper Evident Seals
Product Code Number(s)	5224 BK, 5224 BLACK REMOVE, 5224 GN, 5224 LB, 5224 LV, 5224 OR, 5224 RD, 5224 WT, 5224 YL
Purpose of Product	To show tampering and/or access to critical equipment and supplies.
Range of Applications for Product	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.
Key Specifications of Product	<ul style="list-style-type: none"> • Numbered • Seal and broken by hand • Cannot be resealed once broken • Variety of colors available

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	100 per box
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Seal crash carts, emergency boxes, and other areas where evidence of tampering is critical
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Make sure aperture is ready to accept seal. 2. Place arrow completely through aperture slot/opening. 3. Using thumb and fore finger, apply pressure and close. <ol style="list-style-type: none"> a. You should hear a click to ensure proper seal engagement. b. The seal has now been installed. 4. To remove, twist in any direction to break the shaft/chevron of the seal.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	Follow the policy of your facility for recording date, staff, and other information critical to management of critical inventory.
Special Warnings and Cautions	If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory.
Disposal	N/A

Reprocessing Instructions	
Point of se	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A

Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Compatible with steam sterilization. Standard cycle: 132 °C (270 °F) for 4 minutes
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com