

| Date | Tester Initials | Sonocheck Result Left Side | Sonocheck Result Center | Sonocheck Result Right Side | TOSI Result | Lumcheck Result | Tank Degassed | Bath Temp. | Program Cycle Tested | Action Comments |
|------|-----------------|----------------------------|-------------------------|-----------------------------|-------------|-----------------|---------------|------------|----------------------|-----------------|
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |

Testing of the ultrasonic is done daily or each day the ultrasonic is used

SonoCheck Pass = P
 (Vial turns from blue/green to yellow)

SonoCheck Fail = F
 (Vial stays blue/green in color)

Size of tank determines number of SonoChecks used for testing

SMALL: UP TO 5 LITERS MEDIUM: 5-20 LITERS LARGE: ABOVE 20 LITERS

