

Instructions for Use: Brush Rack

Brand Name of Product	Brush Rack
Generic Name of Product	Brush Rack
Product Code Number(s)	BR002
Intended Use	Brush rack holds the brushes when not in use.
Range of Applications for Product	To hold the brushes and allow the liquid from the bristles to collect in the drip tray.
Key Specifications of Product	<ul style="list-style-type: none"> • 16W inches x 35L inches • Stainless Steel • Wall-mount • 12 angled threaded pegs • Built-in drip tray

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To hold the brushes and allow the liquid from the brushes to collect in the drip tray.
Preparation for Use	1. Securely mount the brush rack to the wall surface with three screws. (Fig.1) .

<p>Diagrams (drawings, pictures)</p>	<div style="text-align: center;">  <p>Figure 1</p> </div> <div style="display: flex; justify-content: space-around; align-items: center;">   </div> <p style="text-align: center;">Figure 2</p>
---	---

Steps for Use of Product	1. Place the eyelet of the brush on the peg to hang on the rack allowing the bristles of the brush to drip into the drip tray. (Fig. 2).
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Wipe the brush rack and drip tray clean with a neutral pH detergent and warm water.
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com