

## Instructions for Use: Brush Storage Rack & Brush Soak Tray

|  |   |
|--|---|
| <b>Brand Name of Product</b>             | Brush Storage Rack & Brush Soak Tray  |
| <b>Generic Name of Product</b>           | Brush Storage Rack & Brush Soak Tray  |
| <b>Product Code Number(s)</b>            | BR001, BRS03  |
| <b>Intended Use</b>                      | The storage rack holds the brushes to allow the bristles to be submerged in cleaning solution in the soaking tray.  |
| <b>Range of Applications for Product</b> | To hold the brushes and allow bristles to be submerged in the soaking tray.   |
| <b>Key Specifications of Product</b>     | <ul style="list-style-type: none"> <li>• Clear Plexiglas</li> <li>• Two holes for mounting.</li> <li>• BR001- Brush Storage Rack             <ul style="list-style-type: none"> <li>○ 10 slots per rack</li> </ul> </li> <li>• BRS03- Brush Soaking Tray</li> </ul> |

| Shipping & Storage                            |  |
|---|--|
| <b>Shipping Conditions &amp; Requirements</b> |  |
| <b>Storage Conditions</b>                     |  |
| <b>Packaging Contents</b>                     |  |
| <b>Shelf Life</b>                             |  |

| Instructions for Using Product           |  |
|--|--|
| <b>Description of Use(s)</b>             | To hold the brushes and allow bristles to be submerged into the soaking tray.  |
| <b>Preparation for Use</b>               | <ol style="list-style-type: none"> <li>1. Mount the storage rack to the wall surface with two screws.</li> <li>2. Mount the brush soak tray(s) below the storage rack at different heights to accommodate both shorter and longer brushes stored on the same rack.</li> <li>3. Use two screws per soaking tray to be mounted to the wall surface.</li> </ol> |
| <b>Diagrams (drawings, pictures)</b>     |  |
| <b>Steps for Use of Product</b>          | <ol style="list-style-type: none"> <li>1. Place the bristles of the brush end into the soaking tray.</li> <li>2. Place the eyelet of the brush on the hook to hang on the storage rack.</li> </ol>   |
| <b>Interpretation of Test Results</b>    |  |
| <b>Contraindications of Test Results</b> |  |
| <b>Documentation</b>                     |  |
| <b>Special Warnings and Cautions</b>     |  |
| <b>Disposal</b>                          |  |

| Reprocessing Instructions                   |   |
|---|---|
| <b>Point of Use</b>                         |   |
| <b>Preparation for Decontamination</b>      |   |
| <b>Disassembly Instructions</b>             |   |
| <b>Cleaning – Manual</b>                    | <ul style="list-style-type: none"> <li>• Remove the soaking tray from the wall and empty the cleaning solution.</li> <li>• Clean the storage rack and soaking tray with neutral pH detergent and warm water.</li> </ul> |
| <b>Cleaning – Automated</b>                 |   |
| <b>Disinfection</b>                         |   |
| <b>Drying</b>                               |   |
| <b>Maintenance, Inspection, and Testing</b> |   |
| <b>Reassembly Instructions</b>              |   |
| <b>Packaging</b>                            |   |
| <b>Sterilization</b>                        |   |
| <b>Storage</b>                              |   |
| <b>Additional Information</b>               |   |
| <b>Related Healthmark Products</b>          |   |

|  |  |
|--|--|
| <b>Other Product Support Documents</b> | ProSys™ Brochure, ProSys™ Price List   |
| <b>Reference Documents</b>             |  |
| <b>Customer Service Contact</b>        | Healthmark Industries Company, Inc.<br>18600 Malyn Blvd.<br>Fraser, MI 48026<br>1-586-774-7600<br><a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a><br><a href="http://hmark.com">hmark.com</a> |