

Instructions for Use: Disposable Bulb & Bladder

Brand Name of Product	Disposable Bulb & Bladder
Generic Name of Product	Disposable bulb & bladder for sigmoidoscopes.
Product Code Number(s)	D30200-210, SBB-01-5
Purpose of Product	For insufflation during rigid sigmoidoscopy. Single use helps eliminate risk of cross-contamination between patients. Cost effective compared to lengthy reprocessing time required for proper cleaning of reusable bulb & bladder.
Range of Applications for Product	Rigid sigmoidoscopy
Key Specifications of Product	<ul style="list-style-type: none"> • Single use. • 100% Latex Free. • No filter required. • Available packaged as non-sterile or sterile.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	Four (4) years from date of manufacture.

Instructions for Using Product	
Description of Use(s)	For insufflation during rigid sigmoidoscopy procedure.
Preparation for Use	N/A
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	Slide black neoprene tubing from insufflation bulb over connection adapter.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	Single use only. Do not reprocess or sterilize.
Disposal	Treat used bulb & bladder as a biohazard and dispose of in accordance with facility policy of contaminated items.

Reprocessing Instructions	
Point of use	N/A
Preparation for decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, inspection, and testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

