



**Instructions for Use: Disposable Bulb & Bladder**

<b>Brand Name of Product</b>	Disposable Bulb & Bladder
<b>Generic Name of Product</b>	Disposable bulb & bladder for sigmoidoscopes
<b>Product Code Number(s)</b>	D30200-210, D30200-210S
<b>Purpose of Product</b>	For insufflation during rigid sigmoidoscopy. Single-use helps eliminate risk of cross contamination between patients. Cost effective compared to lengthy reprocessing time required for proper cleaning of reusable bulb & bladder.
<b>Range of Applications for Product</b>	Rigid sigmoidoscopy
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Single-use. 100% Latex Free.</li> <li>• No filter required.</li> <li>• Available packaged as non-sterile or sterile.</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	4 years from date of manufacture

**Instructions for Using Product**

<b>Description of Use(s)</b>	For insufflation during rigid sigmoidoscopy procedure.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	N/A
<b>Steps for Use of Product</b>	Slide black neoprene tubing from insufflation bulb over connection adapter.
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	Single-use only. Do not reprocess or sterilize.
<b>Disposal</b>	Treat used bulb & bladder as a biohazard and dispose of in accordance with facility policy of contaminated items.

**Reprocessing Instructions**

<b>Point of use</b>	N/A
<b>Preparation for decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, inspection, and testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com